

evo+ Quick Start Guide

Welcome to evo+! Your Steeper engineer has now installed the evo+ app on your device and your profile has been configured based on your requirements for home control. Follow these 3 steps to get started.

1

To start evo+ press this icon on your device screen.



2


As it begins to load, you'll see this homepage. You don't need to press anything else here. This may take a few seconds.



3

When evo+ has fully loaded, you'll be shown your individual configuration complete with your chosen colours and favourite icons.

Now you're ready to go!



evo+ Troubleshooter

We hope you don't encounter any problems with your evo+, however if you do, please first try out these quick troubleshooting tips.

1

Check all switches have working batteries, all controllers are turned on, all connections are secure and that you are connected to the right network.

2

Check the device e.g. TV is turned on, and that it is working from its original remote control.

3

Close and restart the evo+ app.

4

Turn off and restart your iOS device, then restart the evo+ app.

If these steps don't solve the problem you are experiencing, please email us on assist@steepergroup.com or call **0113 207 0449** to speak to our friendly team who will be more than happy to help.