



- FEATURING:**
- IntelliWeave™ composite technology
  - Water friendly EnviroShell™
  - Integrated pyramid

## WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty Inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Composites and Adapters
- CPI Sock
- Foot Shell

**WARRANTY INSPECTION SCHEDULE FOR BREEZE: ANNUALLY.**

## Technical Assistance / Emergency Service 24-7-365

College Park's regular office hours are Monday thru Friday, 8:30 AM – 5:30 PM (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.



**MADE IN THE USA** CE  
383 INS BE TIS 140815

**COLLEGE PARK INDUSTRIES**  
27955 College Park Dr. Warren, MI 48088 USA

**EC REP EMERGO EUROPE**  
Molenstraat 15, 2513 BH, The Hague, Netherlands


©2014 College Park Industries, Inc. All rights reserved. College Park is a registered trademark, with Breeze, IntelliWeave, EnviroShell and FootHorn as trademarks of College Park Industries, Inc.

PACKAGE CONTENTS	TOOLS RECOMMENDED
(1) Breeze™ Foot (1) Foot Shell (1) CPI Sock	4mm Hex Key FootHorn™



## FOOT COMPONENTS

This diagram is to help familiarize you with the unique parts of the Breeze. These parts are referenced in the instructions and used when speaking with a technical service representative.




### COMPONENTS

- 1 Integrated Pyramid
- 2 Housing (Ankle)
- 3 Foreheel with IntelliWeave™ technology
- 4 CPI Sock (not shown)

## DRAIN HOLE

A thin area is located on the bottom of the EnviroShell™. Create a draining path for water by doffing the footshell and drilling or cutting a 1/4" (6mm) diameter hole in the designated location.



**BOTTOM OF SHELL**

## GENERAL INSTRUCTIONS

### MOUNTING

Use only high quality endoskeletal components.

### ASSEMBLY AND DISASSEMBLY (for sock replacement)

Use the FootHorn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

### STATIC ALIGNMENT

For optimal function have the patient's weight balanced between the forefoot and heel. The load line divides the foot at the 1/3 heel lever and 2/3 toe lever. The Breeze was designed with a 3/8" heel rise. Adjustments will need to be made when used without shoes.

### DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE
<b>Firmer Toe Response</b>	Plantarflex the Breeze or move load line posterior
<b>Softer Toe Response</b>	Dorsiflex the Breeze or move load line anterior
<b>Firmer Heel Response</b>	Dorsiflex the Breeze or move load line anterior
<b>Softer Heel Response</b>	Plantarflex the Breeze or move load line posterior

### EXPOSURE TO MOISTURE

To maintain the lifespan of the EnviroShell and CPI sock, always dry foot thoroughly after exposure to moisture. Rinse and dry if exposed to salt water or chlorinated water.

### COSMETIC FINISHING

College Park does not recommend cosmetically finishing the prosthesis if it will be exposed to moisture.

### WARNING

- Under certain circumstances, use in wet environments may increase the risk of slip and fall.
- Remove all shoes and sock and allow to dry in a cool, dry location.
- Do not expose this product to corrosive materials or pH extremes.
- Contaminants such as dirt and the use of lubricants or powder, may effect the function of the CPI Sock and lead to noise.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.