

TECHNICAL INSTRUCTIONS



- FEATURING:**
- IntelliWeave™ composite technology
 - An active heel for natural motion
 - 300 lb. weight limit

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups according to the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- Composites and Adapters
- CPI Sock
- Foot Shell

WARRANTY INSPECTION SCHEDULE FOR CELSUS: ANNUALLY.

Technical Assistance / Emergency Service 24-7-365

College Park's regular office hours are Monday thru Friday, 8:30 AM – 5:30 PM (EST). After hours, an emergency Technical Service number is available to contact a College Park representative.



MADE IN THE USA
126 INS CS TIS 141215

COLLEGE PARK INDUSTRIES
27955 College Park Dr. Warren, MI 48088 USA

EMERGO EUROPE
Molenstraat 15, 2513 BH, The Hague, Netherlands


©2014 College Park Industries, Inc. All rights reserved. College Park and Gait Matching are registered trademarks, with Celsus and FootHorn as trademarks of College Park Industries, Inc.

PACKAGE CONTENTS	TOOLS RECOMMENDED
(1) Celsus Foot (1) Foot Shell (1) CPI Sock	4mm Hex Key FootHorn™



FOOT COMPONENTS

This diagram is to help familiarize you with the unique parts of the Celsus. These parts are referenced in the instructions and used when speaking with technical service representative.



COMPONENTS

- 1

Integrated Pyramid
- 2

Housing (Ankle)
- 3

Fulcrum Pad
- 4

Foreheel
- 5

Heel Wedge (permanent)
- 6

Heel Spring

GENERAL INSTRUCTIONS

MOUNTING

Use only high quality endoskeletal components.

ASSEMBLY AND DISASSEMBLY (for sock replacement)

Use the FootHorn to don and doff the foot shell. Remove the CPI Sock and replace as needed. Any further disassembly or modification of components will void the warranty.

STATIC ALIGNMENT

For optimal function have the patient’s weight balanced between the forefoot and heel. The load line divides the foot at the 1/3 heel lever and 2/3 toe lever.

DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE
Firmer Toe Response	Plantarflex the Celsus or move load line posterior
Softer Toe Response	Dorsiflex the Celsus or move load line anterior
Firmer Heel Response	Dorsiflex the Celsus or move load line anterior
Softer Heel Response	Plantarflex the Celsus or move load line posterior

GAIT MATCHING® GUIDELINES

Contact College Park Technical Service if you have questions about modifying the gait match.

WARNING

- Do not expose this product to corrosive materials, salt water or pH extremes.
- Contaminants such as dirt and the use of lubricants or powder, may effect the function of the CPI Sock and lead to noise.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.