COLLEGE PARK  
**VICTOR™**

*pneumatic knee*

*joelho pneumático · pneumatisches Kniegelenk · rodilla neumática*  
*ركبة هوائية · genou pneumatique · joelho pneumático*  
*пневматическое колено · 气动膝盖*

**college park**

TECHNOLOGY *for the* HUMAN RACE

*technical instructions*


*instruções técnicas · technische anleitung*  
*instrucciones técnicas · التعليمات الفنية*  
*Instructions techniques · Instruções técnicas*  
*Технические инструкции · 技术说明书*

# TABLE OF CONTENTS

<b>EN</b>	Technical Instructions .....	<b>1 - 5</b>
<b>BR</b>	Instruções Técnicas .....	<b>6 - 10</b>
<b>DE</b>	Technische Anleitung .....	<b>11 - 15</b>
<b>ES</b>	Instrucciones Técnicas .....	<b>16 - 20</b>
<b>EU</b>	التعليمات الفنية .....	<b>21 - 25</b>
<b>FR</b>	Instructions Techniques .....	<b>26 - 30</b>
<b>PT</b>	Instruções Técnicas .....	<b>31 - 35</b>
<b>RU</b>	Технические инструкции .....	<b>36 - 40</b>
<b>ZH</b>	技术说明书.....	<b>41 - 45</b>

COLLEGE PARK  
**VICTOR™**

## KEY COMPONENTS


1. Pyramid
2. Upper Assembly
3. Lower Assembly
4. 30mm Pylon Receiver
5. 4mm Clamp Screw (Torque 10 N-m)

## PACKAGE CONTENTS

- (1) Victor Knee
- (1) 4mm Hex Key

## CLEARANCE

1	Overall height	8.0 in (20.4 cm)
2	Dome to knee center	1.1 in (2.9 cm)
3	Dome to tube end contact	6.3 in (16.0 cm)

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

# ADJUSTMENTS

All adjustments can be made using a 4mm Allen wrench.


## STANCE PHASE

### Load Dependence

Controls the brake sensitivity. When load is applied, the knee will not bend until the load is displaced. Turning the load adjustment clockwise will increase the amount of load required to initiate braking.

Factory setting = Screw is at minimum (most brake sensitivity).

*Note: The factory setting for load dependence should generally not require adjustment. However, it may be necessary in cases such as patients with higher body weights, more aggressive gait, or alignment differences.*

### Stance Adjustment


Controls the amount of stance flexion. Adjustments can be made to alter the amount of stance flexion during loading response.

Factory setting = Turn screw clockwise to the point where resistance is felt, then back out about 1/8 turn (counterclockwise).

Readjustment after continued use – The screw should generally not require adjustment; however, play may develop over time. Turn screw 1/8 clockwise to reduce play.

*Note: Stance adjustment works together with load dependence to control the brake function. If decreasing the stance adjustment, it may be necessary to also increase the load dependence adjustment.*

**⚠ Caution:** Do not over tighten. May lead to brake sticking which can cause unstable gait.

## SWING PHASE

### Flexion Resistance

Controls the rate of flexion.

### Extension Resistance

Controls the rate of extension.

Factory setting = both screws at minimum (counterclockwise).

Begin by adjusting the flexion resistance or heel rise first, then the extension resistance to control terminal impact.

*Note: After adjusting the swing phase resistance, it may be necessary to readjust the brake.*

**⚠ Caution:** Flexion and extension must be possible at all settings.

## DYNAMIC ADJUSTMENTS

Note: The effect of a dynamic adjustment can be felt with as little as 1/8 - 1/4 turn.

### STANCE PHASE

#### *Load Dependence*

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Brake function too strong; Swing too difficult to initiate	Increase Load Dependence (brake less sensitive)	Turn <b>L</b> clockwise	
Brake function not sufficient; Knee too unstable during stance	Decrease Load Dependence (brake more sensitive)	Turn <b>L</b> counterclockwise	

#### *Stance Adjustment*

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Not enough flexion	Increase Stance Flexion	Turn <b>S</b> counterclockwise	
Too much flexion	Decrease Stance Flexion	Turn <b>S</b> clockwise	

### SWING PHASE

#### *Flexion Resistance*

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Flexion too fast or Excessive heel rise	Increase Resistance (close valve)	Turn <b>F</b> clockwise	
Flexion too slow or Heel rise not sufficient	Decrease Resistance (open valve)	Turn <b>F</b> counterclockwise	

#### *Extension Resistance*

SYMPTOM	DESIRED RESULT	SCREW ADJUSTMENT	
Extension too fast or Impact too abrupt	Increase Resistance (close valve)	Turn <b>E</b> clockwise	
Extension too slow or Not sufficient	Decrease Resistance (open valve)	Turn <b>E</b> counterclockwise	

## WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

*High patient weight or activity level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.*

- Knee Assembly, Extension Stop Bumpers

**Warranty inspection schedule for College Park Victor:** Six months, then annually.

### WARNING

- Flexion and extension must be possible at all settings.
- Avoid pinching hazards! Do not place fingers near the flexing area of the knee.
- Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to: noise, sudden loss of function, brake release sticking, etc.
- Do not disassemble the knee. Contact College Park to arrange a repair or replacement.
- Do not expose this product to moisture, i.e. fresh water, corrosive materials, salt water or pH extremes. If the knee encounters moisture, wipe it dry using a lint-free cloth.
- Contaminants such as dirt and the use of lubricants or powder may affect the function of the knee brake and lead to failure.
- Do not use compressed air to clean the knee as it can push dirt inside of the knee.
- Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

## TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST).


After hours, an emergency Technical Service number is available to contact a College Park representative.

### LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

### CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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©2019 College Park Industries, Inc. All rights reserved.