

First 3 Liner

The First 3 liner is a clear, soft silicone interface liner with integrated thicker distal walls to provide additional cushioning. It is suitable for below knee (transtibial) amputees up to mobility class 2. It features the unique self-gliding outer surface to allow for easier fitting and removal without the need for alcohol sprays or talcum powder. The liner has a 2-zone matrix which reduces shear forces on the underlying soft tissue. As there is no fabric cover there is no restriction of the liner's radial elasticity.


CES OSS LO1: Read the First liner clinical evaluation summary on the Steeper website.




Features and Benefits



Residual Limb Suitability



Even residual limb shape, decreasing slightly distally, good soft tissue cover




Residual limb bonier, rather poor soft tissue cover and/or with severe scars

More conical residual limb shape, good soft tissue cover

Ordering Information

Size	5+3cm Matrix Part Number	10+3cm Matrix Part Number
18cm	NMF94318	NMF94118
20cm	NMF94320	NMF94120
22cm	NMF94322	NMF94122
23.5cm	NMF94323	NMF94123
25cm	NMF94325	NMF94125
26.5cm	NMF94326	NMF94126
28cm	NMF94328	NMF94128
30cm	NMF94330	NMF94130
32cm	NMF94332	NMF94132
34cm	NMF94334	NMF94134
36cm	NMF94336	NMF94136
38cm	NMF94338	NMF94138
40cm	NMF94340	NMF94140
42cm	NMF94342	NMF94142
45cm	NMF94345	NMF94145

Size Selection





- Determining the size**
1. Measure the circumference of the residual limb 4cm above the distal end of the residual limb.
 2. Attention: Bony structures (e.g. end of the tibia) must be included in the circumference. If the distal tibia is not located within the above measurement, measure higher than 4cm above the distal end of the residual limb.
 3. The circumference minus 0-1 liner size determines the size to be ordered.

- Determining the matrix**
1. Measure the distance from the distal end of the residual limb to the head of the fibula.
 2. Attention: The matrix must not go above the head of the fibula.

- Notes on trimming**
1. The matrix must not be cut.
 2. The trim line for a lower leg prosthesis should lie at least 10cm above the region of the socket edge.
 3. Use the Össur pinking scissors for trimming (Part No. N1950006).

Technical Information

Material	Silicone, Transparent
Matrix	Proprioceptive 2-zone matrix, internal 5+3cm and 10+3cm
Length	Approx 35cm
Amputation Level	Transtibial
Distal Connector	Pin Adapter
Items Included	Silicone glue, Össur washing and drying stand
Warranty	6 months



Replaceable With:



Össur Liner Protect 3C
Page 144

Össur Liner Relax 3C
Page 154