SILITOUCH

1S711 / 1S712

1S721 / 1S722

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EN Instruction for use

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NO Bruerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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	Silitouch liners – 1S711 – 1S712 – 1S721 – 1S722 <i>Instructions for use</i> Read before use	1S71199 Rev. 3 2024-08
---	--	---

Train the patient in how to use the device correctly, and explain to them the advice from § 2, 10, 12, 14, 15, 16, 17, 18, 19. Give the patient these instructions when delivering the device.

1. INCLUDED ITEMS

Part Description*	Part Number* model - thickness - size
Silitouch - LINER TT Aloe locking 3MM - xx	1S711-3-xxx
Silitouch - LINER TT Aloe locking 6MM - xx	1S711-6-xxx
Silitouch - LINER TT Classic locking 3MM - xx	1S712-3-xxx
Silitouch - LINER TT Classic locking 6MM - xx	1S712-6-xxx
Silitouch - LINER TT Aloe w/o locking 3MM - xx	1S721-3-xxx
Silitouch - LINER TT Aloe w/o locking 6MM - xx	1S721-6-xxx
Silitouch - LINER TT Classic w/o locking 3MM - xx	1S722-3-xxx
Silitouch - LINER TT Classic w/o locking 6MM - xx	1S722-6-xxx

* The xx represent the sizes in cm in the part description, and in mm in the part number.

2. DESCRIPTION

The device is a silicone liner which acts as an interface between an amputee's residual limb and the socket of their prosthesis. This device is intended for multiple use by a single patient.

3. PROPERTIES

Length	400 to 450 mm depending on size
---------------	---------------------------------

Model	Distal thickness (mm)	Thickness (mm)	Available sizes (cm)	Distal attachment	Integrated matrix
Silitouch - LINER TT Aloe locking 1S711	17.5 to 18.5	3	16 – 18 – 20 – 22 – 23.5 – 25 – 26.5 28 – 30 – 32 – 34 – 36 – 40 – 45	Yes	10 cm
	19 to 20	6	16 – 18 – 20 – 22 – 23.5 – 25 – 26.5 28 – 30 – 32 – 34 – 36 – 40 – 45	Yes	10 cm
Silitouch - LINER TT Classic locking 1S712	17.5 to 18.5	3	16 – 18 – 20 – 22 – 23.5 – 25 – 26.5 28 – 30 – 32 – 34 – 36 – 40 – 45	Yes	10 cm
	19 to 20	6	20 – 22 – 23.5 – 25 – 26.5 – 28 30 – 32	Yes	10 cm
Silitouch - LINER TT Aloe w/o locking 1S721	14.5 to 16	3	16 – 18 – 20 – 22 – 23.5 – 25 26.5 – 28 – 30 – 32 – 34 – 36 – 40	No	10 cm
	16 to 17.5	6	16 – 18 – 20 – 22 – 23.5 – 25 26.5 – 28 – 30 – 32 – 34 – 36	No	10 cm
Silitouch - LINER TT Classic w/o locking 1S722	14.5 to 16	3	16 – 18 – 20 – 22 – 23.5 – 25 26.5 – 28 – 30 – 32 – 34 – 36 – 40	No	10 cm
	16 to 17.5	6	16 – 18 – 20 – 22 – 23.5 – 25 26.5 – 28 – 30 – 32 – 34 – 36	No	10 cm

4. MECHANISM OF ACTION

The liners have a pressure distributing and shock absorbing function, for improved comfort and protecting the residual limb against moderate impacts and shearing.

The prosthesis suspension is provided by a clutch lock for liners with a distal attachment system, and by vacuum retention (valve with barrier sheath) for other liners.

5. INTENDED USE/INDICATIONS

Silitouch liners are intended for use on a lower-limb prosthetic device and fitted on the patient by a certified prosthetist (health professional). They are indicated for transtibial amputees with low to high activity levels (1 to 4).

Liners with aloe vera are recommended for patients with more sensitive skin.

6. CLINICAL BENEFITS

- Compression of the residual limb (preventing bruising, reducing fluctuations in volume of the residual limb)
- Comfort, flexibility
- Low risk of irritation and sensitisation
- Adhesion to the residual limb
- Firm hold on the prosthesis

7. ACCESSORIES AND COMPATIBILITY

Liners with a distal attachment have a standard distal attachment system, with a diameter M10 thread.

⚠ Use a lock compatible with the liner's distal attachment system.

⚠ Copolymer gel cushions/pads are not compatible with these liners.

Part Number	Part Description	Use
1S109-80	80 mm silicone distal cushion	Silicone cushion for volume compensation and/or protecting fragile areas on the residual limb
1S109-100	100 mm silicone distal cushion	
1S109-120	120 mm silicone distal cushion	
1X300	STEPLINE hood	A textile hood for volume variation compensation, which is compatible with liners without a distal attachment system (available in sizes 20 to 60 cm)
1X301	STEPLINE hood for DA	A textile hood for volume variation compensation, which is compatible with distal attachment liners (available in sizes 25 to 50 cm)


8. CHOOSING THE SIZE

1. Measure the circumference of the residual limb 4 cm from the distal end.
2. Choose the liner size corresponding to the closest size below this measurement.

⚠ Do not choose a liner bigger than the measurement taken.

⚠ Never go down more than one liner size. This could cause excessive compression.

⚠ The matrix for distal attachment liners should end under the fibula head, and not protrude past it. This could impede knee flexion.

9. TRIMMING (BY A CERTIFIED PROSTHETIST ONLY)

⚠ Do not trim under the matrix.


⚠ Excessive trimming could reduce the adhesion between the liner and the limb, affecting the fastening of the prosthesis.

⚠ Use a pair of scissors in good condition.

⚠ In case of fraying after trimming, stick the fabric back on using a biocompatible silicone glue.

10. FITTING THE LINER

1. Check the condition of the liner, and make sure that there are no solid particles in the liner, which could cause a skin reaction. Before first using the liner, clean the inside.
2. Cover open wounds to prevent direct contact between the wound and the liner.
3. Turn the liner over and place it on the end of the residual limb, preventing air bubbles from forming.
4. For distal attachment liners, check that the **locking system is centred and aligned** in the longitudinal axis of the residual limb.
5. Carefully unroll the liner onto the residual limb.


⚠ Avoid creases in the area behind the knee.

⚠ Avoid pulling on the upper end of the liner. That could cause increased tension on the edge of the liner and the underlying skin and cause a skin reaction.

6. If some creases appeared when putting on the liner, repeat the procedure.

⚠ Take care not to damage the silicone with your fingernails or any other sharp object when handling the liner.

11. TROUBLESHOOTING (FOR THE CERTIFIED PROSTHETIST)

Problem	Solution
Breakdown of fabric or silicone in the anterior-distal area of the tibia.	Check the socket fit.
Fabric wear along the proximal end of the socket.	- Smooth the edge of the socket. - Check that the liner is not rubbing on the socket during flexion/extension. - Add an edge lining, such as leather, onto the edge if necessary.
Liner piston movements in the socket.	Check the socket fit.
Residual limb and liner rotating inside the socket.	Check the socket fit.
Silicone detaching from the fabric on the edge of the liner.	Replace the liner.
Upper edge of the liner rolling up the leg.	Check that the patient is wearing a liner the right size.
Excessive sweating	Check the socket fit, and if necessary, choose a liner of a different size thickness or style.
Skin irritation all over the residual limb.	- Review the cleaning procedures with the patient. - Check whether the patient has changed cleaning product. - Check the socket fit, and if necessary, choose a liner of a different size thickness or style. - Have a physician check that it is not an allergic reaction (patch test).
Skin irritation along the proximal edge of the liner.	- Re-explain the instructions for putting on the liner: it must be unrolled, and not pulled. - Trim the edge of the liner in an undulating shape.

12. WARNINGS

⚠ If the packaging is damaged, check the integrity and cleanliness of the device.

⚠ A faulty product can cause loss of adhesion. In this case, stop using the liner and consult your certified prosthetist.

⚠ Store the liner out of the reach of children. They could choke by putting it over their mouth and nose.

⚠ Make sure to keep the device's packaging bearing the CE label.

13. CONTRAINDICATIONS

- Excessive soft tissue
- Residual limb cannot carry loads
- Poor hygiene
- Allergy/intolerance to material
- Irritated skin, subject to dermatitis

14. SIDE EFFECTS

- Excessive perspiration
- Skin reaction

⚠ In case of a skin reaction, consult your certified prosthetist and physician.

Any serious incident that occurs and that relates to the device must be reported to the manufacturer and to the competent authority in the Member State where the user is based.

15. CARE

Have the fit and condition of the liner (edges, tearing/any irregularities, condition of silicone, etc.) checked by the certified prosthetist at least every 6 months, and replace it if necessary.

16. CLEANING THE LINER AND SKIN

⚠ Regularly clean the inside and outside of the liner, and disinfect it at least once a week.

Cleaning

- Machine: put the liner in a laundry bag, and wash at 30°C with a mild detergent at slow spin 
- By hand: turn the liner inside-out (fabric on the inside), and wash in warm water (30°C) with a balanced pH soap, turn the liner the right way round (fabric on the outside), and then wipe down the inside and outside with a towel

Leave the liner to dry in the open air. Do not dry it on a heat source. **Do not leave inside-out or twisted during drying.** Do not put the liner back on when still wet.

⚠ Do not use deodorising products, dyes, softeners, stain removers, whitening agents or bleaches, or any other non-recommended chemical. This could damage the liner.

⚠ Do not handle the liner close to unprocessed residue of fibreglass, carbon fibre or other particles which might not be eliminated by normal cleaning. The presence of these residues on the silicone could cause a skin reaction.

Skin care

Wash the residual limb daily with a pH neutral soap, free from deodorising products and dyes.

17. ENVIRONMENTAL CONDITIONS

Store and use the liner at ambient temperature. Do not store near a heat source.

18. DISPOSAL

The liner comprises silicone (with or without aloe vera) and a synthetic fibre fabric. The device and packaging must be disposed of in accordance with current local or national regulations.

19. SERVICE LIFE

The liner undergoes natural wear from continuous use. Its service life depends on its use. It is recommended to renew your liner at least every 6 months.

20. DESCRIPTION OF THE SYMBOLS

	Manufacturer		Identified risk		CE marking and year of 1 st declaration		Single patient, multiple use
---	--------------	---	-----------------	---	--	---	------------------------------

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PROTEOR SAS