


# GEWA One By Abilia

**GEWA One is a remote control that helps users manage appliances in their own homes.**

The GEWA One can manage multiple devices including; TV, music system, lamps, recliner chairs, profiling beds, door, curtain and window openers. By inserting a SIM card, the GEWA One also works as a mobile phone enabling the user to keep in touch with friends, family and quickly reach a carer through either a call or text message. For ease of use, the user can choose whether to control the device through the touch screen or navigate with an external switch, and the specially designed symbols, colours and voice recordings can assist the user in quickly finding the right button.



## EASY TO CONTROL


Users only need to tap the screen lightly to select a button and activate an appliance, and the GEWA One can utilise larger buttons for a greater target area. Alternatively, if the user cannot touch the screen they can use a variety of external switches, where each button is highlighted in turn on the screen until the one needed is recognised, upon which the user presses the switch again to select the button. For easier navigation, symbols, different colours and recorded speech on the buttons can be easily programmed and the number of buttons that a page displays can be customised.

## DESIGNED FOR A RANGE OF USERS

The GEWA One is suitable for a wide range of users including; those who have difficulties with operating a conventional remote control due to motor/cognitive/visual disabilities. Those who have difficulties with operating a conventional smartphone or don't want to use their smartphone for environmental control. Those who need to be able to use their remote control sitting in a wheelchair or lying in bed and those who don't like to be dependent on people in their environment.

## SIMPLE INSTALLATION

Our team of experienced engineers can quickly install the GEWA One into the home with the minimum of fuss. The device is easy to customise either on the device or online as the entire process is completed on the myAbilia web portal. Within myAbilia there is a backup copy of the users GEWA One setup making some support functions accessible remotely and quickly without the need for an engineer to visit.



## TECHNICAL INFORMATION

### MOBILE UNIT

- Screen size - 6.2"
- Size - 166 x 83 x 13,3 mm
- Weight - 267 g
- OS - Android 9

### MULTIBOX

- Size - 108 x 67 x 22 mm
- Weight - 153 g
- Connections: 3.5 mm, Micro-USB
- Control signals: IR, Z-Wave Radio

### INCLUDED WITH THE DEVICE;

- Mobile unit, restricted to GEWA One
- Charger mobile unit including USB-cable
- GEWA Multibox
- Charger GEWA Multibox incl USB-cable
- Headphones, 3,5 mm incl USB-C adapter
- USB-cable USB-C to USB-A
- SIM-tray installation tool
- User manual
- Voucher with license key

## CONTACT INFORMATION

- Please contact Steeper Group via phone on 0113 207 0449 or e-mail [assist@steepergroup.com](mailto:assist@steepergroup.com) for more information, to request a demo and to place an order.



Steeper Group,  
Unit 3 Stourton Link,  
Intermezzo Drive,  
Leeds,  
LS10 1DF

Tel: 0113 207 0449  
Email: [assist@steepergroup.com](mailto:assist@steepergroup.com)

[www.steepergroup.com](http://www.steepergroup.com)

STPAT032