

Smartbox Grid Pad Pro Specifications



The powerful communication aid with environmental control, flexible mounting and innovative sound.

- Three sizes to suit your needs
- Built in switch ports and high quality speakers
- Radio and infrared environment control

Powered by Grid 3, our award winning AAC software with symbol and text communication, computer control, interactive learning and more.

Features	10s	12	15
Operating system	Windows 10 Pro	Windows 10 Pro	Windows 10 Pro
Software	Grid 3	Grid 3	Grid 3
Processor	Intel® Core™ i5-7Y54, 1.2GHz Base, 3.2GHz Turbo	Intel® Core™ i5-7Y54 Processor, up to 3.20 GHz	Intel® Core™ i5-7Y54 Processor (1.2GHz, 3.2GHz Burst)
Memory	8GB	8GB DDR3	8GB
Hard drive	256GB SSD	SSD SATA III, 128GB	SSD SATA III, 256GB
Video card	Intel HD Graphics	Intel HD Graphics	Intel HD Graphics
Dimensions	260x189x60mm	355x223x45mm 14x8.8x1.8inch	408x259x42mm 21x10x1.7"
Weight	1.45kg 3.2lbs	2.7kg 5.9lbs	3kg 6.6lbs
Battery life	Up to 7 hours	Up to 15 hours (10 hours with eye gaze)	Up to 14 hours (9 hours with eye gaze)
Tablet camera	2MP front, 5MP rear	2MP front, 5MP rear	2MP front, 5MP rear
Screen	10.1" Full HD display	12.5" Full HD display	15.6" Full HD display and 4.6" rear-facing second screen
Access	Touchscreen, switch, USB input	Touchscreen, switch, USB input	Touchscreen, switch, USB input
Connectivity	Wi-fi, Bluetooth	Wi-fi, Bluetooth	Wi-fi, Bluetooth
Ports	2xUSB 2.0, SD Card Slot	2xUSB 2.0, SD Card Slot	2xUSB 2.0, SD Card Slot
Buttons	Power, Volume	Power, Volume	Power, Volume
Switch access	2x3.5mm switch ports	2x3.5mm switch ports	2x3.5mm switch ports
Mounting plate	Rehadapt and Daessy compatible	Rehadapt and Daessy compatible	Rehadapt and Daessy compatible
Environment control	GEWA infared, Z-Wave radio/ EasyWave	GEWA infared, Z-Wave radio/ EasyWave	GEWA infared, Z-Wave radio/ EasyWave
Wheelchair mounting	REHADapt and Daessy plate with table stand	REHADapt and Daessy plate with table stand	REHADapt and Daessy plate with table stand
Warranty	3 years SmartCare (UK only)	3 years SmartCare (UK only)	3 years SmartCare (UK only)

Smartbox Grid Pad Pro Specifications



The powerful communication aid with environmental control, flexible mounting and innovative sound.

- Three sizes to suit your needs
- Built in switch ports and high quality speakers
- Radio and infrared environment control

Powered by Grid 3, our award winning AAC software with symbol and text communication, computer control, interactive learning and more.

Features	10s	12	15
Processor	Intel® Core™ i5-7Y54, 1.2GHz Base, 3.2GHz Turbo	Intel® Core™ i5-7Y54 Processor, up to 3.20 GHz	Intel® Core™ i5-7Y54 Processor (1.2GHz, 3.2GHz Burst)
Hard drive	256GB	128GB	256GB
Dimensions	260x189x60mm 10x7.4x2"	355x223x45mm 14x8.8x1.8"	408x259x42mm 21x10x1.7"
Weight	1.45kg 3.2lbs	2.7kg 5.9lbs	3kg 6.6lbs
Battery life	Up to 7 hours	Up to 15 hours (10 hours with eye gaze)	Up to 14 hours (9 hours with eye gaze)
Tablet camera	2MP front	2MP front	2MP front, 5MP rear
Screen	10.1" Full HD display	12.5" Full HD display	15.6" Full HD display and 4.6" rear-facing second screen

	10s	12	15
Warranty (3 years)*	✓	✓	✓
Second screen			✓
Tablet camera	✓	✓	✓
WiFi	✓	✓	✓
Bluetooth	✓	✓	✓
Grid 3 software	✓	✓	✓
Antimicrobial screen	✓	✓	✓
Mountable	✓	✓	✓
Screen	✓	✓	✓
USB Access	✓	✓	✓
Windows 10 Pro	✓	✓	✓
8GB Memory	✓	✓	✓
Eyegaze	✓	✓	✓
Switch access	✓	✓	✓

*Smart care warranty is only available in the UK