



Assistive Technology
ACTIV 500

 **steeper**
Creating life's turning points, together

ACTIV 500

Independence at the touch of a button

The ACTIV 500 is a simple to operate handset, designed to enable independent living. The system can carry out numerous simultaneous functions for comprehensive home automation, including intercoms, door entry, phones, TV and DVD.



Home Control

Facilitating a large range of everyday tasks the ACTIV 500 can be configured to operate over 200 functions around the house, granting newfound freedom to users, whatever their disability.



Flexibility

Featuring a large, easy press button and two connections for other input switches, as well as speech prompts. The ACTIV 500 has a flexible configuration, and easy to use set up making it an ideal choice for a wide variety of clients.



Radio and Infrared Controls

The ACTIV 500 is compatible with existing Steeper radio control and mainstream infrared peripherals for easy upgrading.



Key Features

- Simple to use scanning controller with large, back-lit screen.
- Configurable to the user's needs, including auditory prompting, scanning tones and icons.
- Multi-lingual with recorded speech.
- Can be used as a basic communication device with selectable pre-recorded voice phrases.
- Powerful infrared learning technology can control anything requiring a remote control.