

Universal Palm Protector



Features and Benefits

- Fleece-lined, lightweight palm cover
- Designed to allow full finger movement whilst reducing finger contractures and skin damage/breakdown in the palmar region
- Easily laundered
- Easily applied using hook and loop fastening straps

Indications for Use

- Post surgery/injury positioning
- Hypertonic contractures
- Rheumatoid arthritis

Ordering Information

Description	Size	Product Code
Right	Universal	PPR
Left	Universal	PPL

CMC Brace



Features and Benefits

- Support for the thumb CMC joint in the treatment of osteoarthritis and medial subluxation. The shaped strap wraps around the thumb CMC joint and the first metacarpal for targeted CMC joint support
- Made from Neoprene

Indications for Use

- CMC joint arthritis
- CMC instability
- Tendinitis
- Pain associated with repetitive motion

Ordering Information

Size	Palm Circ. (cm)	Product Code	
		Left	Right
Small	16.5 - 17.5	CMC/SL	CMC/SR
Medium	17.5 - 20	CMC/ML	CMC/MR
Large	20 - 24	CMC/LL	CMC/LR

Thumb Spica - Neoprene



Features and Benefits

- Neoprene thumb spica featuring a removable flexible stay which restricts movement at the CMC and MCP joints, leaving the IP joint free
- Warmth and compression
- Removable flexible stay

Indications for Use

- CMC joint arthritis
- CMC instability
- Tendinitis
- Pain associated with repetitive motion
- Osteoarthritis
- Painful/unstable 1st MCP joint
- Rheumatoid arthritis
- Thumb strains
- Sprains and post-cast removal

Ordering Information

Size	Wrist Circ. (cm)	Product Code	
		Left	Right
Small	14 - 15	TS/BE/S/L	TS/BE/S/R
Medium	15 - 16.5	TS/BE/M/L	TS/BE/M/R
Large	16.5 - 18	TS/BE/L/L	TS/BE/L/R
X-Large	18 - 19	TS/BE/XL/L	TS/BE/XL/R