

Collamed® Knee Brace




Features and Benefits

- Long soft brace with the possibility of extension and flexion limitation
- The compressive spacer-knit can improve proprioception
- Easy to handle thanks to wrap design
- Silicone dot coating supports a secure fit on the leg
- Medi Airtex+ material with intelligent thermoregulation provides a high level of comfort
- Special knit technology at the hollow of the knee allows ease of flexion
- The physioglide hinge mimics the roll-glide movement of the knee and supports a secure fit on the leg
- Comprehensive set of wedges for flexible extension and flexion limitation including 5° hyperextension wedge. The hinge design enables the wedges to be changed without using a tool

Indications for Use

- Injuries to the collateral ligament (post traumatic/postoperative)
- Injuries to the cruciate ligament (post traumatic/postoperative)
- Instability of the knee
- Meniscus injuries
- Osteoarthritis of the knee



Note: Measure 15cm above the centre of the patella.

Ordering Information

Size	Thigh Circ. (cm)	Product Code
X-Small	31 - 37	824/1
Small	37 - 45	824/2
Medium	45 - 51	824/3
Large	51 - 57	824/4
X-Large	57 - 65	824/5
XX-Large	65 - 74	824/6
XXX-Large	74 - 84	824/7

Rebound® Post-Op Knee Brace



Features and Benefits

- Simple, tool-free design
- Multi-purpose Smart-Fit™ packaging
- OneTouch™ easy to adjust ROM hinge with lock-out capability
- Colour-coding and numbered frame for patient education
- Telescoping frame from 47cm - 68cm (Adult) and 39.5cm - 54cm (paediatric)
- Quick-fit, gripping buckles for easy application
- The brace is lightweight and the edges are protected to improve patient comfort

Indications for Use

- ACL, PCL, MCL, LCL repairs/injuries
- Osteochondral repairs
- Meniscal repairs
- Condylar fractures
- Tibial plateau fractures
- Sprains/strains of the knee
- Patellar tendon repairs
- HTOs (High Tibial Osteotomy)

Note: Measure 15cm above the centre of the patella.

Ordering Information

Size	Description	Thigh Circ. (cm)	Product Code
Universal	Cool	41 - 81	B-8381480
Universal	Full Foam	41 - 81	B-8382480