

## Wrist and Thumb Spica - Neoprene



### Features and Benefits

- Restricted wrist and thumb motion without limiting finger motion
- Easy fit and comfortable
- Lightweight and padded
- Contours to the anatomy of the wrist and thumb
- Breathable Neoprene material
- Hook and loop fastening

### Indications for Use

- Gamekeeper's thumb
- Post-cast healing
- Soft tissue injuries
- Selective scaphoid injuries



### Ordering Information

Size	Wrist Circ. (cm)	Product Code	
		Left	Right
Small	15 - 18	WTS/S/L	WTS/S/R
Medium	18 - 20	WTS/M/L	WTS/M/R
Large	20 - 23	WTS/L/L	WTS/L/R
X-Large	23 - 25	WTS/XL/L	WTS/XL/R

## Comfyrene® Contour Hand Orthosis



### Features and Benefits

- The Comfyrene® Contour Hand Orthosis has unique side flaps that hold the wrist in place allowing for full hand contouring
- C-bar and thumb support with included strap
- 4 straps come as standard
- Unique butterfly designed side straps hold wrist and forearm in place
- Contours the whole hand, wrist and forearm
- Instead of removable terrycloth covers, the Comfyrene splint inserts are covered with soft neoprene and headliner material, which creates a slim profile for everyday comfort
- Splints are adjustable and re-adjustable using the Comfy™ patented reinforced spine
- Straps may be positioned anywhere on the splint for optimal fit
- Contracture management

### Indications for Use

- Use after sports injuries
- Removal of casts
- Sports medicine
- Arthritic patient's support of the weakened extremity

### Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code	
			Left	Right
Adult	17.8+	27.9	HCO-101-CP-LB-LT	HCO-101-CP-LB-RT

## Comfy Air Hand Orthosis



### Features and Benefits

- The Comfy Air Hand Orthosis is an exceptionally useful splint for the severely contracted hand
- It has an adjustable air bladder that allows gradual extension of the MP, DIP and PIP joints, as well as an adjustable wrist for customised fit
- Ingeniously designed by a therapist, one simply inserts the unfitted splint with the bladder under the fingers and then pumps air to inflate the bladder and increase finger extension. (Headliner material only)
- Ambidextrous

### Indications for Use

- Ideal for the severely contracted hand
- Inflatable air bladder allows for gradual MP, DIP and PIP joint extension
- Adjustable wrist
- Provides low pressure stretch therapy
- Hand air bladder can be used separately without full orthosis as needed
- Can fit into even the most tightly clenched hand

### Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	HA101-A

## Comfy Hand-Wrist Finger Orthosis



### Features and Benefits

- Provides support and positioning for the weak or deformed wrist, hand or fingers
- Serves as an excellent resting splint to prevent trauma to joints and positions to increase ROM
- Side wings prevent ulnar and radial deviation
- Arm cuffs easily adjust to width of arm
- Terrycloth
- Headliner
- Ambidextrous

### Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Carpal tunnel syndrome

**Note:** Paediatric sizes also available, see **Page 216.**

### Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult Small	15.2 - 17.8	27.9	H101-AS
Adult Large	20.3+	34.3	H101-AL