

## ComfySplints™ Rest Hand



### Features and Benefits

- Supports the weak or deformed hand with ultra-comfort
- Removable finger separator included
- Super soft and removable navy headliner cover is machine washable (laundry bag included)
- Treated with antimicrobial barrier for skin breakdown prevention
- Patented reinforced spine can bend and adjust as needed (no tools or heat required)

### Indications for Use

- Arthritis
- CVA
- Sprain/strains



### Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code	
			Left	Right
Adult Small	13.3 - 15.9	26.7	RH-101-AS-LEFT	RH-101-AS-RIGHT
Adult	15.9 - 18.4	29.2	RH-101-LEFT	RH-101-RIGHT
Adult Large	18.4+	34.3	RH-101-AL-LEFT	RH-101-AL-RIGHT

## ComfySplints™ Large Pan Hand Thumb Orthosis



### Features and Benefits

- Designed to support the thumb without stressing the web space between the index finger and the thumb known as the "Thenar Eminence"
- Soft straps accommodate sensitive skin, oedema and fragile skin
- Removable Terrycloth material covers are fully machine washable
- Fully and easily adjustable at all joints to make a custom fit without the use of tools or heat guns
- Side wings prevent ulnar/radial deviation
- Ambidextrous

### Indications for Use

- Sprains/Strains/Fractures
- Sports injuries
- Neurological impairments
- Contractures
- CP
- Wound care



### Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	LPHT-101