

## ComfySplints™ Spring Loaded Goniometer Hand Orthosis



### Features and Benefits

- Goniometer works in the same manner as serial casting by gradually adjusting the amount of extension
- Features a spring that gently pushes extremity to the selected degree of extension
- Use the goniometer dial to set the range of flexion and extension from 0 to 180° in 10° increments
- Allen key supplied
- Blocks flexion or extension while allowing motion in opposite direction
- Can be set in a desired static position
- Re-adjustable custom fit of wrist, MP, PIP and DIP joints
- Ambidextrous

### Indications for Use

- Contractures of the wrist and hand
- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Stroke

### Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	SGH-101



Push Ortho Thumb  
 Brace CMC  
**Page 58**