

Diaped Duosoft Therapeutic



Features and Benefits

- Constructed with a multi-density material combination of plastazote and cellular urethane, these advanced insoles offer superb resilience, cushioning and shock absorption
- The insoles are designed to fit into both traditional footwear and 'additional depth' orthopaedic shoes
- The addition of an increased metatarsal support in the anatomical design is designed to facilitate off-loading of the metatarsal head area

Indications for Use

- For individuals with diabetes, arthritis, or other conditions that result in sensitive feet
- Metatarsalgia
- Pre-ulcerated sites
- Plantar fasciitis
- Foot deformities
- Morton's neuralgia
- Forefoot accommodation
- Skin breakdown

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
2 - 4.5	34 - 37	SB2335
5 - 6.5	36 - 38	SB2336
7 - 8.5	39 - 42	SB2337
9 - 10.5	43 - 45	SB2338
11 - 13	46 - 48	SB2339

Diaped Duosoft Plus



Features and Benefits

- Therapeutic insoles for insensate feet
- The dual-layer design of Duosoft Plus consists of a resilient, shock absorbing urethane base and a PORON 94® top cover
- The slow rebound quality of the material cradles the foot and allows enhanced redistribution of body weight over the foot's plantar surface
- Soft, comfortable and allows the feet to remain cool, dry and odourless
- A Vygel heel accompanies the dual-layer formation to provide optimal shock attenuation at heel strike

Indications for Use

- For individuals with diabetes, arthritis, or other conditions that result in sensitive feet
- Metatarsalgia
- Pre-ulcerated sites
- Plantar fasciitis
- Foot deformities
- Morton's neuralgia
- Forefoot accommodation
- Skin breakdown

Ordering Information

Size	Shoe Size (UK)	Product Code
A	2 - 4	SB2402
B	5 - 7	SB2404
C	8 - 9	SB2406
D	10 - 12	SB2408
E	13 - 15	SB2410