

Foot-Up®



Foot-Up shown with shoeless wrap



Foot-Up shown with footwear



Features and Benefits

- Foot-Up is a lightweight ankle orthosis that offers dynamic support for drop foot or similar complaints
- The revolutionary orthosis provides visible improvement in gait by supporting the foot the moment it is raised
- Lightweight and difficult to see under clothes
- Constructed from a padded breathable material, it is comfortable to wear with or without footwear (shoeless wrap required for use without footwear)

Indications for Use

- Peroneal neuropathies
- Drop foot
- Dysfunctional gait

Styles Available



Ordering Information

| Size | Ankle Circ. (cm) | Product Code | |
|-----------|------------------|--------------|---------|
| | | Black | Beige |
| 1 Medium | 18 - 21 | 07810-1 | 07811-1 |
| 2 Large | 22 - 26 | 07810-2 | 07811-2 |
| 3 X-Large | 27 - 33 | 07810-3 | 07811-3 |

Shoeless Wrap*

| Size | Foot Circ. (cm) | Product Code | |
|-----------|-----------------|--------------|---------|
| | | Black | Beige |
| 1 Small | <20 | 27810-1 | 27811-1 |
| 2 Medium | 20 - 24 | 27810-2 | 27811-2 |
| 3 Large | 24 - 28 | 27810-3 | 27811-3 |
| 4 X-Large | >28 | 27810-4 | 27811-4 |

Plastic Insert

| Description | Product Code |
|-------------|--------------|
| Black | 17810 |
| Beige | 17811 |

- = Size of Foot-Up (1 - 3)
- = Size of Shoeless Wrap (1 - 4)

Example: To order a black size 1 (medium) Foot-Up with a size 2 (medium) Shoeless Wrap, the code would be **07810-12**.

Foot-Up with Shoeless Wrap Bundle

| Description | Product Code |
|-------------|--------------|
| Black | 07810-■● |
| Beige | 07811-■● |

Note: *Shoeless Wrap comprises of foot wrap only and must be ordered with complete Foot-Up kit unless the patient already has one.