

## Mini TLC™ Hip Abduction Brace



### Features and Benefits

- Prefabricated
- Material: MPE
- Pelvic band adjustable for hip development and pelvic circumference with quick release buckle closure
- Thigh cuff bivalve design adjustable to thigh proximal and distal measurements
- Breathable, washable liner absorbs perspiration and reduces heat build up
- Universal hip joint, TLC hip joint with infinite abduction/adduction and flexion/extension settings
- Available in 6 sizes

### Indications for Use

- Arthroplasty
- Hip dislocation or potential high risk dislocation
- Hip dysplasia - anterior or posterior
- Hip management/immobilisation
- Hip revision
- Hip surgery
- Postoperative
- Preoperative
- Stabilise, align and reinforce hip



### Ordering Information

Size	Waist Circ. (cm)	Proximal Thigh Circ. (cm)	Distal Thigh Circ. (cm)	Product Code
XX-Small	40.6 - 45.7	29.2 - 34.3	20.3 - 25.4	TLCMINI2S
X-Small	45.7 - 50.8	31.8 - 36.8	22.9 - 27.9	TLCMINIXS
Small	50.8 - 55.9	34.3 - 39.4	25.4 - 30.5	TLCMINISM
Medium	55.9 - 61	36.8 - 41.9	27.9 - 33	TLCMINIMD
Large	61 - 66	39.4 - 44.4	30.5 - 35.6	TLCMINILG
X-Large	66 - 71.1	41.9 - 47	33 - 38.1	TLCMINIXL

## Jobskin® Leg Gaiters Standard Lining



### Features and Benefits

- External orthosis with structural inserts to immobilise the affected joint and hold the limb in controlled extension, or the best possible position for rehabilitation
- Providing a prolonged stretch, immobilisation and contracture management to help maintain the gains achieved by manipulation, corrective surgery or reconstructive procedures
- For use on one or more limbs and can be used by infants through to adulthood
- Designed to meet the individual's functional needs rather than to treat an individual condition or pathology
- Ready-to-wear

### Indications for Use

- Certain skin conditions e.g. Eczema
- Vascular/Circulation problems
- Oedema
- Weight fluctuations
- Epilepsy
- Contracture management of the lower limb

**Note:** This product is non-returnable. A pattern must be stipulated when ordering. For available patterns, see **page 221**.

### Ordering Information

Size	Top Circ. (cm)	Bottom Circ. (cm)	Length (cm)	Product Code
X-Small	18 - 25	13 - 20	20	LG090220
			25	LG090225
Small	24 - 31	15 - 22	20	LG090320
			25	LG090325
			30	LG090330
Medium	31 - 38	18 - 25	30	LG090430
			35	LG090435
			40	LG090540
Large	39 - 46	21 - 28	45	LG090545
			50	LG090550
			45	LG090645
X-Large	46 - 53	24 - 31	50	LG090650