

Oasis LSO



Features and Benefits

Oasis is our most popular bivalve spinal system and ideal for patients who require maximum support, but cannot tolerate most rigid LSO braces.

- The traction closure system provides easy compression adjustments
- Revolutionary spacer fabric offers increased breathability, heat and moisture removal, less bulk and longer life span than trilaminates
- Flexibility of closing mechanism ensures an excellent fit across various body shapes
- Effective between L2 - L5

Indications for Use

- Burst fracture
- Post-op lumbar spine
- Stabilisation
- Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- Spinal instability
- Post-op disc surgery
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours or inflammation
- Post-op laminectomy

Ordering Information

Size	Waist Circ. (cm)	Hip Circ. (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OALS-XS
Small	76.2 - 91.4	88.9 - 104.1	OALS-SM
Medium	86.4 - 101.6	99 - 114.3	OALS-MD
Large	96.5 - 111.8	109.2 - 124.5	OALS-LG
X-Large	106.7 - 121.9	119.4 - 134.6	OALS-XL
XX-Large	116.8 - 132	129.5 - 144.8	OALS-XXL
XXX-Large	127 - 142.2	139.7 - 154.9	OALS-XXXL

Terrapin LSO Bivalve Brace



Features and Benefits

- Bivalve design allows for easy donning and doffing
- Provides abdominal compression
- Adjustable A/P control
- Removable and washable liner
- Made from high density polyethylene
- Standard option has a front height of 10" (25.4cm)
- Short option has a front height of 8½" (21.6cm)

Indications for Use

- Post operative spinal immobilisation
- Compression fracture
- Spinal instability
- Chronic back pain
- Degenerative, herniated or bulging disc

Ordering Information

Size	Hip Circ. (cm)	Hip Circ. (inch)	Product Code		
			Standard	Short	Pendulous Abdomen
1	76.2 - 86.4	30 - 34	LSO-1	-	-
2	86.4 - 94	34 - 37	LSO-2	-	-
3	94 - 104.1	37 - 41	LSO-3	LSO-3S	-
4	104.1 - 114.3	41 - 45	LSO-4	LSO-4S	LSO-4P
5	114.3 - 124.5	45 - 49	LSO-5	-	LSO-5P