

Formfit® Walker



Features and Benefits

- Low, wide rocker bottom promotes a stable, natural gait
- Shock-absorbing sole for reduced impact at heel strike
- Breathable, seamless liner
- Available in black

Indications for Use

- For foot and ankle conditions that may benefit from motion restriction such as:
 - Grade II or III ankle sprains
 - Stable fractures
 - Soft tissue injuries
 - During postoperative or rehab phase

Ordering Information

Size	Shoe Size (UK)	Product Code
X-Small	1.5 - 3.5	WO200BLK
Small	4 - 6.5	WO400BLK
Medium	7 - 10	WO600BLK
Large	10.5 - 12	WO800BLK
X-Large	12.5 - 14.5	W1000BLK

Rebound® Air Walker



Features and Benefits

- The Rebound Air Walker features linked medial and lateral air bladders with a single inflation point to provide balanced compression. These air bladders are shaped to cradle the Achilles, holding the heel in place for extra stability
- Moisture wicking soft goods and ventilated panels prevent microbial build up and helps to keep the skin dry and comfortable
- Straps fasten in front for easy application, adjustment and removal
- EVA material provides greater shock absorption in the heel
- Fully-integrated air pump and release valve provides easy, tool-free inflation and deflation for fine-tuned, consistent compression
- Flex Edge™ over-moulding conforms to each patient's unique anatomy to help limit the formation of pressure points
- Rocker bottom promotes a natural, stable gait, whilst a semi-soft toe bumper that can be trimmed as required provides added protection

Indications for Use

- Grade 2 and 3 sprains
- Stable fractures
- Postoperative use
- Trauma
- Rehabilitation use

Ordering Information

Size	Men's Shoe Size (UK)	Women's Shoe Size (UK)	Product Code
X-Small	2 - 4	3.5 - 5.5	B-242900001E
Small	4 - 7	5.5 - 8.5	B-242900002E
Medium	7 - 10.5	8.5 - 12	B-242900003E
Large	10.5 - 12.5	12 - 14	B-242900004E
X-Large	12.5 - 15	14 - 16.5	B-242900005E

Rebound® Diabetic Walker



Features and Benefits

- Uniquely designed for effective treatment of plantar foot ulcers, the Rebound Diabetic Walker enables easy access and re-fitting during wound care, whilst improving the patient's mobility during the healing process by including a rocker bottom outsole for smooth roll-over
- Effective tri-laminate insole with removable shock-absorbing hexagons that promote ulcer offloading and a Plastazote top layer that evenly distributes the pressure
- Anatomically designed shell that minimises shear forces by immobilising the foot and ankle at 90° whilst providing protection
- Flex Edge technology for relief of pressure points, along with a patent-pending pressure relief valve to prevent over-inflation
- Positioning portals allow the clinician to check the fit of the walker at the heel bottom, in order to ensure appropriate positioning of the foot in the walker
- Soft and breathable machine washable liner for hygiene promotion and prevention of infection

Indications for Use

- Plantar foot ulcers

Ordering Information

Size	Men's Shoe Size (UK)	Women's Shoe Size (UK)	Product Code
Small	4.5 - 7	5.5 - 8.5	B-242500002
Medium	7 - 10.5	8.5 - 12	B-242500003
Large	10.5 - 12.5	12 - 14	B-242500004
X-Large	12.5 - 15	14 - 16.5	B-242500005

VACOped



Features and Benefits

- Inspired by ski boots, the VACOped features a unique design to enable patients to weight-bear earlier, enabling a faster return to daily activities
- Inner vacuum cushion provides equal and firm support by adjusting itself precisely to the shape of every foot and lower leg and allowing adjustments to be made for swelling
- Adjustable ROM bar allows the boot to be fixed or permit movement dependent on the indication
- Durable, yet lightweight German-engineered plastic shell provides the injury site with the necessary stability and allows an increased freedom of movement in comparison to a plaster cast
- Multiple washable and exchangeable liners promote hygiene

Indications for Use

- Achilles tendon ruptures
- Ankle fractures
- Metatarsal fractures
- Calcaneus fractures
- Ankle replacements
- Arthrodesis

Ordering Information

Size	Shoe Size (UK)	Calf Circ. (cm)	Ankle Circ. (cm)	Product Code
Small	2.5 - 6.5	23 - 39	14 - 29	VPED-01-4-S
Medium	6.5 - 10.5	35 - 51	22 - 40	VPED-01-4-M
Large	10.5 - 14	40 - 58	30 - 45	VPED-01-4-L