

A-Boot Universal



Features and Benefits

- Ambulation version (A-Boot) comes with removable, soft, non-slip sole
- Antimicrobial surface protector to shield against mould and mildew, stains, bacteria, fungus and algae
- Super soft fleece lined headliner cover is removable and machine washable
- A toe extension is available to accommodate all foot sizes
- Rotator bar controls hip and leg rotation
- Boots position heel in floating position to prevent pressure areas through pressure reduction
- Comes with extender kit to accommodate for wider or oedematous feet
- Universal sizing to reduce stock holding

Indications for Use

- Wider or oedematous feet
- Ankle contractures
- Heel ulcerations

Note: Paediatric sizes also available, see **Page 231**.

Ordering Information

Size	Product Code
Universal	A-BOOT

More Information
 Wear and Care Instructions

Universal PRAFO 650 Kodel



Features and Benefits

- Designed for recumbent or ambulant protection of the heel and foot
- Zero pressure or shear force at the heel
- Universal sizing to reduce stock holding
- Resists plantarflexion contractures developing
- Custom modifications are possible
- One version fits both left and right
- Calibrated adjustment of plantar/dorsiflexion
- Integral walking base - designed for ambulation
- Washable and replaceable kodel liner system
- Adjustable foot, toe and calf sections
- Aluminium heel connector for strength
- Rotation bar resists internal/external rotation when used in bed
- Suitable for day and night use

Indications for Use

- Pressure ulcers
- Can be used to resist the development of contracture

Note: Paediatric sizes also available, see **Page 230**.

Ordering Information

Size	Product Code	Accessories	Product Code
Universal	PRFO	Universal liner	PRAFO-LINER

Universal PRAFO APU (Adjustable Posterior Upright)



Features and Benefits

- Designed for recumbent or ambulant protection of the heel and foot
- Zero pressure or shear force at the heel
- Resists plantarflexion contractures developing
- Custom modifications are possible
- One version fits both left and right
- Calibrated adjustment of plantar/dorsiflexion
- Integral walking base - designed for ambulation
- Washable and replaceable kodel liner system
- Adjustable foot, toe and calf sections
- Aluminium heel connector for strength
- Rotation bar resists internal/external rotation when used in bed
- Universal sizing to reduce stock holding

Indications for Use

- Pressure ulcers
- Can be used to resist the development of contracture

Note: Paediatric sizes also available, see **Page 230**.



Ordering Information

Size	Product Code	Accessories	Product Code
Universal	PRFO/APU	Universal liner	PRAFO-LINER

Universal PRAFO Dual Action



Features and Benefits

- The pad and strap version allows bulky wound dressings to be accommodated and can allow superior foot control when the PRAFO® is to be used in walking training and other circumstances when mobility is the primary focus
- Universal sizing to reduce stock holding
- Bi-planar adjustments with a screwdriver
- Opportunity to evaluate and train gait
- Ambulatory or recumbent usage
- All the benefits of the 650 and 650 APU - with increased application possibilities
- Calibrated adjustment of plantar/dorsiflexion
- Integral walking base - designed for ambulation
- Washable and replaceable kodel liner system
- Adjustable foot, toe and calf sections
- Aluminium heel connector for strength
- Rotation bar resists internal/external rotation when used in bed
- Universal liners are also available to order - use code PRAFO-LINER

Indications for Use

- The Dual Action accommodates and helps control varus or valgus conditions of the ankle/foot complex, as well as allowing the practitioner a more precise and simplified method of setting the ankle/foot complex in the sagittal plane

Note: Paediatric sizes also available, see **Page 230**.

Ordering Information

Size	Product Code	Accessories	Product Code
Universal	PRFO/DA	Universal liner	PRAFO-LINER