

PRIM Dorsi Assist

Features and Benefits

- Consists of two separated pieces, a supramalleolar strap over the ankle and a non-slip support union to fix between the tongue and shoelaces by means of a central anchorage hook type. Both parts are joined by an elastic strap, which allows variable amounts of dorsiflexion assist
- On the back of the ankle band, there is a gel pad which is cushioned for increased comfort on the achilles tendon
- Discrete and flexible, it takes up no space in the shoe making it easy to fit
- To be fitted when lacing the shoe
- Has the function to produce dorsiflexion on the take-off phase of gait, indicated in flaccid palsy

Indications for Use

- Drop foot
- Peroneal neuropathies
- Dysfunctional gait



Plantar Band			
Size	Measurement (cm)	Product Code	
		Left	Right
Small	17 - 21	ATX11-SL	ATX11-SR
Medium	21 - 25	ATX11-ML	ATX11-MR
Large	25 - 29	ATX11-LL	ATX11-LR



Ordering Information

Textile Drop Foot Orthosis		
Size	Measurement (cm)	Product Code
X-Small	13 - 17	ATX01-XS
Small	17 - 21	ATX01-S
Medium	21 - 25	ATX01-M
Large	25 - 29	ATX01-L

Accessories	
Description	Product Code
Shoe Support	ATX12
Elastics Small	ATX13-S
Elastics Medium	ATX13-M
Elastics Large	ATX13-L
Metal Hooks (2 Pack)	ATX14