

Kids Feet In Motion



Features and Benefits

- These brightly coloured pairs consist of one blue (left) and one right (red) to help children learn whilst providing a functional solution
- Medium density EVA orthosis specifically designed for children's feet
- A flat base on the standard orthotic allows for modification with wedges and components
- Insoles are available as standard or with a pre-posted heel of 3° or 5°

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Planovalgus feet
- Postural related symptoms

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code		
		Standard	3°	5°
4	19 - 20	SC3304	SC3330	SC3350
5	21 - 22	SC3305	SC3331	SC3351
6.5	23 - 24	SC3306	SC3332	SC3352
8	25 - 26	SC3307	SC3333	SC3353
9 - 10	27 - 28	SC3308	SC3334	SC3354
11	29 - 30	SC3309	SC3335	SC3355
12 - 13	31 - 32	SC3310	SC3336	SC3356
1 - 2	33 - 34	SC3311	SC3337	SC3357

EZIKID (Peddy) Insoles



Features and Benefits

- Dynamic insole for children
- Modular sandwich-like construction and flexible thermoformable wedge
- Soft sole to support the natural movement of the foot

Indications for Use

- Back pain
- Posture related symptoms
- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Planovalgus feet

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
4	19 - 20	ZKR1-0640-19/20
5	21 - 22	ZKR1-0640-21/22
6.5	23 - 24	ZKR1-0640-23/24
8	25 - 26	ZKR1-0640-25/26
9 - 10	27 - 28	ZKR1-0640-27/28
11	29 - 30	ZKR1-0640-29/30
12 - 13	31 - 32	ZKR1-0640-31/32
1 - 2	33 - 34	ZKR1-0640-33/34
2 - 3	35 - 36	ZKR1-0640-35/36
4 - 5	37 - 38	ZKR1-0640-37/38
6 - 7	39 - 40	ZKR1-0640-39/40

Podo Tech Junior ¾ Length



Features and Benefits

- Designed especially for children
- Easily ground and heat mouldable

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
4	19 - 20	SC3220
5	21 - 22	SC3222
6.5	23 - 24	SC3224
8	25 - 26	SC3226
9 - 10	27 - 28	SC3228
11	29 - 30	SC3230
12 - 13	31 - 32	SC3232
1 - 2	33 - 34	SC3234
2 - 3	35 - 36	SC3236
4 - 5	37 - 38	SC3238
6 - 7	39 - 40	SC3240
8 - 9	41 - 42	SC3242

PediPod Full Length Insoles



Features and Benefits

PediPod insoles are uniquely designed to help children develop healthy foot and body posture as they grow. The development of poor posture while growing potentially disturbs normal walking, running and standing which can lead to recurrent back, knee and foot pain as an adult - often resulting in a lifetime of treatment.

- Anatomical support features in insoles built for children's feet
- Provides comfort and pain relief by deep heel cup and paediatric orthotic insole specific arch support
- The brushed nylon surface is comfortable and allows easy passage of the foot into the shoe
- A full-length insole that can be quickly and easily trimmed to fit most comfort children shoes

Indications for Use

- Aching arches on exercise
- Children's heel pain during and after exercise (Sever's Disease)
- Ankle, calf or leg aches on exercise
- Children's kneecap pain during or after sport or on stairs

Ordering Information

Size	Shoe Size (UK)	Product Code
C	12 - 1	PED03
D	10 - 11.5	PED04
E	8 - 9.5	PED05
F	6 - 7.5	PED06

Peapod Junior



Features and Benefits

- 1.5mm thermoplastic base offers excellent fitting compatibility in most footwear
- Heat mouldable and grindable

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
5	21 - 22	SC2502
6.5	23 - 24	SC2504
8	25 - 26	SC2506
9 - 10	27 - 28	SC2508
11	29 - 30	SC2510
12 - 13	31 - 32	SC2512
1 - 2	33 - 34	SC2514
2 - 3	35 - 36	SC2516

Kraemer Heel Cups



Features and Benefits

- Children's multi-layered prefabricated foot orthosis covered in Multiform A30 shore
- Thermo-mouldable reinforced base and wedge

Indications for Use

- Ligament weakness
- Sprains/strains
- Planovalgus feet

Ordering Information

Shoe Size (UK)	Shoe Size (EU)	Product Code
4	19 - 20	ZKR1-0607-19/20
5	21 - 22	ZKR1-0607-21/22
6.5	23 - 24	ZKR1-0607-23/24
8	25 - 26	ZKR1-0607-25/26
9 - 10	27 - 28	ZKR1-0607-27/28
11	29 - 30	ZKR1-0607-29/30
12 - 13	31 - 32	ZKR1-0607-31/32
1 - 2	33 - 34	ZKR1-0607-33/34
2 - 3	35 - 36	ZKR1-0607-35/36
4 - 5	37 - 38	ZKR1-0607-37/38
6 - 7	39 - 40	ZKR1-0607-39/40