

Universal PRAFO Paediatric



Features and Benefits

- The PRAFO® (Pressure Relief Ankle Foot Orthosis) is a totally adjustable, custom-fitted AFO that can help manage many ankle/foot anomalies, offering a range of adjustable features
- Custom-contoured aluminium heel connector bar helps to control dorsi-plantar flexion by allowing infinite, measurable adjustments
- Ergonomically designed foot plate allows the liner to be changed easily without tools
- Integral, securely fastened vibram walking base provides a shoe-type surface for ambulatory patients
- The design protects the patient around the vulnerable areas of the foot and heel and will remain so during ambulation
- Calf and foot segments can be custom-modified to accommodate patient physiology
- Anti-rotation bar made of malleable aluminium adjusts to control unwanted extremity rotation
- Suitable for day and night
- Universal sizing to reduce stock holding

Indications for Use

- Pressure ulcers
- Can be used to resist the development of contracture

Note: For larger sizing see **Page 131, 132 & 133.**

Ordering Information

Size	Product Code
Universal	PRFO-550

Paediatric A-Boot



Features and Benefits

- Washable, soft fleece-lined cover wicks moisture away from the skin to maintain skin integrity and allows heel to float, thus eliminating pressure of friction on the heel and enhancing blood circulation
- Rubber sole with non-slip surface provides a walking base so patients can walk comfortably
- The splint frame resists plantar flexion while providing a dynamic push toe extension for the correction of foot drop, foot ankle contractures and deformity
- The adjustable toe post relieves pressure on the toes and adjusts to accommodate a variation of foot sizes
- Rotation bar can be positioned to the side to control hip and leg rotation
- An extender kit is available for patients with increased calf and foot circumference and or oedema
- Universal sizing to reduce stock holding

Indications for Use

- For the patient who requires positioning but can also transfer and ambulate

Note: For larger sizing see **Page 130.**

Ordering Information

Size	Foot Length (cm)	Calf Circ. (cm)	Heel to Top of Back (cm)	Product Code
Universal	12.7 - 20.3	<17.9	17.8	A-BOOT-P