

## PeroSupport.tec



### Features and Benefits

- Limits the plantar flexion in the upper ankle joint
- Magnetic closure simplifies opening and closing with one hand
- Great energy storage and return by leveraging carbon stiffness
- Can be trimmed to size with sturdy scissors
- Highly cosmetic

### Indications for Use

- Peroneal paresis
- Weakness of dorsiflexion of the foot
- Drop foot
- Hemiplegia



### Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Height (cm)	Product Code	
				Left	Right
X-Small	2 - 3.5	34 - 36	32.5	197R4/L0	197R4/R0
Small	3.5 - 5.5	36 - 39	33.5	197R4/L1	197R4/R1
Medium	5.5 - 8	39 - 42	34.5	197R4/L2	197R4/R2
Large	8 - 10.5	42 - 45	35.5	197R4/L3	197R4/R3
X-Large	10.5 - 13	45 - 48	37.5	197R4/L4	197R4/R4

Spare Accessories	Size	Product Code
3x pads, 1x magnetic fastener and 10x hook and loop dots	Universal	197R10

## PeroSupport.pro



### Features and Benefits

- Dynamic re-adjustment of the foot
- Restriction of plantar flexion and support of dorsal flexion
- Prevents lateral ankle twisting (supination position)
- Prevents side-slips (lateralisation)
- Great energy storage and return by leveraging carbon stiffness
- Can be trimmed to size with sturdy scissors

### Indications for Use

- Peroneal paresis
- Weakness of dorsiflexion of the foot
- Malalignment of the ankle joint as a result of weakness of dorsiflexion of the foot
- Hemiplegia



### Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Height (cm)	Product Code	
				Left	Right
Small	3.5 - 5.5	36 - 39	34	197R6/L1	197R6/R1
Medium	5.5 - 8	39 - 42	35	197R6/L2	197R6/R2
Large	8 - 10.5	42 - 45	35	197R6/L3	197R6/R3
X-Large	10.5 - 13	45 - 48	36	197R6/L4	197R6/R4

Spare Accessories	Shoe Size (UK)	Shoe Size (EU)	Product Code
Calf pad with fastener, 12 hook and loop dots	3.5 - 8	36 - 42	197R7/1
	8 - 13	42 - 48	197R7/2