

Push Aequi Junior Ankle Brace



Features and Benefits

- Designed for children between the ages of 6 and 12 years old (with shoe size approx. 12UK/30EU - 3.5UK/36EU)
- Stabilises the ankle joint during inversion as well as eversion and neutral and flexed positions
- The freedom of dorsal and plantar flexion is sufficiently maintained to allow normal walking movement
- Optimal functioning of the brace is obtained when it is worn with a sturdy lace-up shoe

Indications for Use

- Acute ankle ligament lesions (both inversion and eversion trauma)
- Stabilisation of conservatively or operatively treated ankle fractures
- (Residual) instability
- Secondary prevention of ankle ligament lesions

Note: For larger sizing see **Page 99**.



Ordering Information

Size	Product Code
Junior	PUS100B

UFO™



Features and Benefits

Introducing an improved prefabricated night splint for the effective treatment of plantar fasciitis and achilles tendonitis. The UFO™ comfortably supports the affected foot in a controlled amount of dorsiflexion, applying a Low Load Prolonged Stretch (LLPS) to the soft tissues of the plantar fascia and achilles tendon. LLPS promotes gentle stretching, healing and remodelling of these tissues, effectively reducing the stress and pain caused by inflammation and tissue shortening. Such support can help break the cycle of injury, pain and re-injury due to otherwise unlimited contraction and deformity of the soft tissues at night.

- The simple variable dorsiflexion design allows the patient to quickly and easily adjust the amount of applied stretch to provide the most effective pain relief
- Adjustable bilateral dorsiflexion assist straps
- CoolFoam™ replaceable liner wicks moisture away from the skin for enhanced hygiene

Indications for Use

- Plantar fasciitis
- Achilles tendonitis
- Contracture management of the foot and ankle

Ordering Information

Size	Footplate Length (cm)	Brace Height (cm)	Product Code
Small	14	16.5	B/3614
Medium	17	25	B/3615
Large	19	26.5	B/3616