

Formfit® Shoulder Brace



Features and Benefits

- Streamlined fit and sporty look
- Breathable and comfortable
- Provides abduction of the shoulder (when the abduction pillow is attached), or external rotation of the shoulder (when the external rotation pillow is attached)
- Available with a customisable aeroplane kit for patients who require precise positioning following shoulder surgery
- Optional axilla pillow available for purchase to help distribute weight over a larger surface area and decrease pressure on the ulnar nerve

Indications for Use

- Immobilisation for rotator cuff repairs
- Bankart repairs
- Capsular shifts
- Glenohumeral dislocations/subluxations
- Soft tissue repairs/strains
- Shoulder instabilities/injuries
- Post-surgical rehabilitation



Ordering Information

Size	Length (cm)	Product Code		
		Regular	External Rotation	With Abduction
Small	Up to 29	SSS-SO	SSS-ER	SSS
Medium	29 - 40	SSR-SO	SSR-ER	SSR
Large	40 - 44.5	-	-	SSL

Neuro-Lux Shoulder II® (Sporlastic)



Features and Benefits

- Freedom of movement allows virtually unrestricted use of affected arm
- Orthosis can be taken on and off with one hand
- Different coloured plastic buckles make it easy to put on correctly
- Improved repositioning of joint (humerus) whilst correcting subluxation and inner rotation
- Stabilisation of shoulder joint
- Prevents secondary injuries of crown cap, tendons, muscles and shoulder nerves
- Improvement of posture through movement of arm pendulum
- Avoidance of shoulder-hand syndrome

Indications for Use

- Flaccid and painful movement of the shoulder and arm musculature after hemiplegia
- Plexus injuries
- Craniocerebral trauma
- Cervical spine trauma
- Peripheral nerve injuries

Ordering Information

Size	Lower Arm Circ (cm)	Chest Circ (cm)	Product Code	
			Left	Right
1	21 - 25	75 - 90	07261/1L	07261/1R
2	21 - 25	90 - 105	07261/2L	07261/2R
3	21 - 25	105 - 125	07261/3L	07261/3R
4	25 - 29	75 - 90	07261/4L	07261/4R
5	25 - 29	90 - 105	07261/5L	07261/5R
6	25 - 29	105 - 125	07261/6L	07261/6R
7	25 - 35	75 - 90	07261/7L	07261/7R
8	25 - 35	90 - 105	07261/8L	07261/8R
9	25 - 35	105 - 125	07261/9L	07261/9R