

## SoftPro® Functional Resting WHFO



### Features and Benefits

- “Bend-to-fit” endoskeleton frame can be hand shaped in the sagittal plane as well as medially, laterally, and rotationally to accommodate most angles of the wrist/hand/and fingers without the use of tools
- The endoskeleton frame is a closed cell foam, covered by a laminated foam/cloth cover that can be laundered as necessary
- Extra long wrist straps hold the wrist in the proper position and applies a gentle stretch
- Maintains and supports the wrist, hand and fingers in a comfortable functional position

### Indications for Use

- Support and protection of a weak or injured wrist, hand, and fingers, including thumb

### Ordering Information

Size	Width of MP Joints (cm)	Length (cm)	Product Code	
			Left	Right
Small	6.5 - 7.5	23	SP31FFR-S-L	SP31FFR-S-R
Medium	7.5 - 9	25	SP31FFR-M-L	SP31FFR-M-R
Large	9 - 10	28	SP31FFR-L-L	SP31FFR-L-R

## ComfySplints™ Deviation Hand Thumb Orthosis



### Features and Benefits

- Unique, below the wrist, swivel that allows for lateral and medial positions to accommodate for ulnar and radial deviation
- Progressively change wrist and finger positions and amount of wrist deviation by adjusting the internal frame and locking the frame's swivel with included hex wrench
- Correct or control wrist deviation motion even when changes in the position of the wrist and hand are made
- Ambidextrous

### Indications for Use

- Arthritic deformities
- Neuromuscular deformities
- Flaccid, weak or spastic hand conditions

### Ordering Information

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	DHT-101