

Standard Knee Brace



Features and Benefits

- The front opening design allows easy application for the arthritic or elderly
- Metal side hinges provide medial/lateral stability while allowing flexion
- Opposing fastenings for upper and lower touch and close straps counter rotation and allow correct fitting
- Removable hinges
- Universal left/right application
- Made from Airprene

Indications for Use

- Medio lateral stability
- Osteoarthritis
- Rheumatoid arthritis

Note: When measuring the knee, ensure it is fully extended. For thigh circumference, measure 15cm above the centre of the patella. For calf circumference, measure 15cm below the centre of the patella.

Ordering Information

Size	Thigh Circ. (cm)	Knee Circ. (cm)	Calf Circ. (cm)	Product Code
X-Small	29 - 34	18 - 22	26 - 31	SKB/CX/BE/XS
Small	34 - 39	22 - 27	31 - 37	SKB/CX/BE/S
Medium	39 - 47	27 - 29	37 - 42	SKB/CX/BE/M
Large	47 - 53	29 - 34	42 - 46	SKB/CX/BE/L
X-Large	53 - 59	34 - 40	46 - 53	SKB/CX/BE/XL




More Information
 Wear and Care Instructions

Stabimed® Knee Brace



Features and Benefits

- Short soft brace with the possibility of extension and flexion limitation
- Easy to handle thanks to wrap design
- Silicone dot coating supports a secure fit on the leg
- Medi Airtex+ material with intelligent thermoregulation provides a high level of comfort
- Special knit technology at the hollow of the knee allows ease of flexion
- The physioglide hinge mimics the roll-glide movement of the knee and supports a secure fit on the leg
- The compressive spacer-knit can improve proprioception
- Comprehensive set of wedges for flexible extension and flexion limitation including 5° hyperextension wedge. The hinge design enables the wedges to be changed without using a tool



Note: Measure 15cm above the centre of the patella.

Indications for Use

- Injuries of the collateral ligaments
- Instability of the knee
- Meniscus injuries
- Osteoarthritis of the knee

Ordering Information

Size	Thigh Circ. (cm)	Product Code
X-Small	31 - 37	826/1
Small	37 - 45	826/2
Medium	45 - 51	826/3
Large	51 - 57	826/4
X-Large	57 - 65	826/5
XX-Large	65 - 74	826/6