



## Swedish Type Lightweight AFO

### Features and Benefits

- Heat mouldable splint providing dorsiflexion assistance and lateral stability for ankle/foot area
- Easily trimmed and low profile with open heel to fit into a shoe easily
- Wide calf strap with touch/close fastening
- Removable/washable padding for comfort and hygiene

### Indications for Use

- Flaccid foot drop after neurological impairment

### Ordering Information

| Size    | Heel - Calf Height (cm) | Footplate Length (cm) | Product Code |           |
|---------|-------------------------|-----------------------|--------------|-----------|
|         |                         |                       | Left         | Right     |
| Small   | 24                      | 26                    | AFO-1S-L     | AFO-1S-R  |
| Medium  | 29                      | 27.5                  | AFO-2M-L     | AFO-2M-R  |
| Large   | 33                      | 28                    | AFO-3L-L     | AFO-3L-R  |
| X-Large | 38                      | 30                    | AFO-4XL-L    | AFO-4XL-R |



## Swedish Type Lightweight AFO - Carbon Reinforced

### Features and Benefits


- Heat mouldable splint providing dorsiflexion assistance
- Easily trimmed with low profile and open heel to fit into a shoe easily
- Wide calf strap with touch and close fastening and removable/washable padding for comfort and hygiene


### Indications for Use

- Flaccid foot drop after neurological impairment

### Ordering Information

| Size    | Heel - Calf Height (cm) | Footplate Length (cm) | Product Code |            |
|---------|-------------------------|-----------------------|--------------|------------|
|         |                         |                       | Left         | Right      |
| Small   | 24                      | 26                    | AFO-1S-LC    | AFO-1S-RC  |
| Medium  | 29                      | 27.5                  | AFO-2M-LC    | AFO-2M-RC  |
| Large   | 33                      | 28                    | AFO-3L-LC    | AFO-3L-RC  |
| X-Large | 38                      | 30                    | AFO-4XL-LC   | AFO-4XL-RC |

 **More Information**  
 Wear and Care Instructions

 **More Information**  
 Patient Information