

# VENUM Universal TLSO PPX

## Features and Benefits

The VENUM Universal TLSO PPX is the ultimate one-size-fits-all semi-rigid comfort brace for the thoraco lumbar-sacral region. This brace reduces pain, whilst providing support and immobilisation to the thoraco lumbar-sacral spine. Featuring unique anterior pads and a highly effective cincher system, the VENUM Universal TLSO PPX is designed to produce maximal intra-abdominal pressure to distract the vertebral bodies of the spine. As the brace is universal in nature, it benefits from fitting every patient type and can be stored easily where limited storage is available.

- Reduced stock holding
- Cincher technology
- Cost effective
- Step-down, streamline design
- Motion restriction of the thoracic and lumbar regions, T7 - L5
- Mechanically advanced single or dual compression closure system provides a 6:1 pull ratio
- Posterior panel is available in independent or central support panel designs which creates a natural support system for the spine that can be customised to provide the perfect fit
- Ergonomic grips can be customised to patient size by extending or retracting cable length
- Made from durable, breathable, soft, hypo-allergenic materials that wick away perspiration
- Fits waist circumferences from 60cm - 157.5cm; optional extender panel adds 25cm to circumference

Universal



**More Information**  
Patient Fitting Guide

## Indications for Use

- Chronic back pain
- Degenerative, herniated or bulging disc
- Spondylolysis
- Spondylolisthesis
- Spinal stenosis
- Facet syndrome
- Tumours
- Spinal instability
- Burst fracture
- Pre or post-op spinal thoraco-lumbar stabilisation
- Pre or post-op disc surgery

**Note:** To fit the VENUM Universal TLSO PPX correctly, please choose the desired size on the slide panel from the charts below.

## Ordering Information

Description	Size	Product Code
VENUM Universal TLSO PPX	Universal	VIPHYBTLSO4PPX
Extender Panel	Universal	STEXT

## VENUM Universal TLSO and Hybrid TLSO Size Chart

Size	Waist Measurement (cm)	Size	Waist Measurement (cm)
X-Small	60 - 65	X-Small/Small	65 - 70
Small	70 - 75	Small/Medium	75 - 80
Medium	80 - 85	Medium/Large	85 - 90
Large	90 - 95	Large/X-Large	95 - 100
X-Large	100 - 105	X-Large/2X-Large	105 - 110
2X-Large	110 - 115	2X-Large/3X-Large	115 - 120
3X-Large	120 - 125	3X-Large/4X-Large	125 - 130
4X-Large	130 - 135	4X-Large/5X-Large	135 - 140
5X-Large	140 - 145	Add extender panel	157.5 - 182.5
Panel fully open	142.5 - 157.5	-	-