

VENUM Universal Soft Hip Abduction Brace



Features and Benefits

The one-size-fits-all VENUM Soft Hip Abduction brace incorporates the most advanced technology and materials. A revolutionary soft hip management system for patients of all shapes and sizes providing maximum comfort and support. The LSO part of the brace operates using a cincher system to maximise control. The belt thigh cuff and hinge adjust to fit the dimension of most patients. The brace applies a dynamic abduction force, so that the femoral head remains in the socket.

- Universal design
- Cincher technology
- Size indicator displays belt size to ensure correct and comfortable fitting
- Joint with adjustable abduction/adduction, flexion/extension
- Single or double-pull cable pulley system
- Stock reduction, saves money and space
- Materials are comfortable, breathable, soft, wick-away perspiration and are durable and resilient
- Waist circumference range: 60cm - 157.5cm
- Extender panel adds an additional 25cm (sold separately)
- Ergonomic arthritic grips can be customised to patient size

Indications for Use

- For patients at risk of hip dislocation, suitable for all levels of instability

Material: 30% Polyurethane 70% Polyamide

Ordering Information

Description	Size	Patient Height (cm)	Product Code
VENUM Soft Hip Abduction Brace	Universal	<188	VHHIPABD
Long Joint*	Universal	>188	TLCHIPJNTLG
Extender Panel	Universal	-	STEXT

***Note:** The Venum Universal Soft Hip Aduction Brace comes with a standard joint for patients of height below 188cm, for patients over 188cm a long joint will be required.

To fit the VENUM Universal Hip Brace correctly, please choose the desired size on the slide panel from the chart below.

VENUM Universal Braces Size Chart

Size	Waist Measurement (cm)	Size	Waist Measurement (cm)
X-Small	60 - 65	X-Small/Small	65 - 70
Small	70 - 75	Small/Medium	75 - 80
Medium	80 - 85	Medium/Large	85 - 90
Large	90 - 95	Large/X-Large	95 - 100
X-Large	100 - 105	X-Large/2X-Large	105 - 110
2X-Large	110 - 115	2X-Large/3X-Large	115 - 120
3X-Large	120 - 125	3X-Large/4X-Large	125 - 130
4X-Large	130 - 135	4X-Large/5X-Large	135 - 140
5X-Large	140 - 145	Add extender panel	157.5 - 182.5
Panel fully open	142.5 - 157.5	-	-