

## Universal Wrist Brace



Universal



### Features and Benefits

- Breathable felt and neoprene mix material
- Aluminium palmer splint and stays can be shaped to fit the wrist with additional hook and loop strap for enhanced fit
- Designed for easy on and off application
- Circumferential seam for improved fit around waist

### Indications for Use

- Arthritis
- Sprain
- Swelling
- Carpal Tunnel

**Note:** Paediatric sizes also available, see **Page 214.**

### Ordering Information

Description	Size	Product Code
Left	Adult	SWB/L/UNIVERSAL
Right	Adult	SWB/R/UNIVERSAL

## Wrist and Thumb Spica



### Features and Benefits

- Two-piece dual compression with malleable aluminium thumb splint to hold thumb in abduction
- 5cm wide elasticated strap
- Made from Neoprene

### Indications for Use

- Osteoarthritis
- Painful/unstable 1st MCP joint
- Rheumatoid arthritis

### Ordering Information

Size	Wrist Circ. (cm)	Product Code	
		Left	Right
Small	12.5 - 15	WTB/CX/BE/S/L	WTB/CX/BE/S/R
Medium	15 - 17.5	WTB/CX/BE/M/L	WTB/CX/BE/M/R
Large	17.5 - 20	WTB/CX/BE/L/L	WTB/CX/BE/L/R
X-Large	20 - 22.5	WTB/CX/BE/XL/L	WTB/CX/BE/XL/R