

Calzamedi Footwear



Features and Benefits

- Extended stiffeners on both medial & lateral sides provide improved stability of the foot
- Lightweight rubber soles contribute to a more energy efficient gait
- Hook and loop fastening or lace styles help children develop independent dressing skills
- Low opening style makes footwear easy to take on and off
- Extra padding around heel stiffeners to improve comfort
- Fully breathable lining reduces skin irritation
- Adaptable sole provides enhanced function
- Leather uppers are durable and breathable
- Padded collars ensure a comfortable fit
- Supplied with removable anatomically shaped footbeds

Indications for Use

- Flexible flat feet
- Hypermobility
- Hindfoot instability
- Sensitive feet
- Children who need extra support during gait

Note: Please see footwear size chart on the next page (page 244-245).



Pair



Black Leather



On Demand



Blue Nubuck



On Demand



Pink Nubuck

Ordering Information

Black Leather			
Fastening Type		Product Code	
Lace		■	● -4200-0101-
Hook and Loop		■	● -4201-0101-
Pink Nubuck		Blue Nubuck	
Fastening Type	Product Code	Fastening Type	Product Code
Lace	■ -4202-0203-	Lace	■ -4202-0202-
Hook and Loop	■ -4203-0203-	Hook and Loop	■ -4203-0202-

■ = Shoe Size (EU)

● = Width Number (4, 6, 8, Anti Varus, or DAFO)

See footwear sizing chart on following pages for more sizing information.

Example: To order a pair of black leather lace up shoes in size 37 with a width of 7.7cm (4 width), the code would be **37-4200-0101-4**.



Calzamedi Footwear Size Chart

Shoe Size (UK)	Shoe Size (EU)	Inner Shoe Length (cm)	Foot Length (cm)	Width 4		Width 6	
				Width Circ. (cm)	Joint Circ. (cm)	Width Circ. (cm)	Joint Circ. (cm)
1	16	10.3	9.2 - 9.7	12.7	4.8	13.7	5.2
1.5	17	11	9.8 - 10.3	13.1	5	14.1	5.4
2.5	18	11.7	10.4 - 10.9	13.5	5.1	14.5	5.5
3	19	12.3	11 - 11.5	13.9	5.2	14.9	5.6
4	20	13	11.6 - 12.1	14.3	5.3	15.3	5.7
5	21	13.7	12.2 - 12.8	14.7	5.4	15.7	5.8
5.5	22	14.3	12.9 - 13.3	15.1	5.5	16.1	5.9
6	23	15	13.4 - 14	15.5	5.6	16.5	6
7	24	15.7	14.1 - 14.7	15.9	5.7	16.9	6.1
8	25	16.3	14.8 - 15.2	16.3	5.8	17.3	6.2
8.5	26	17	15.3 - 15.9	16.7	6	17.7	6.4
9.5	27	17.7	16 - 16.6	17.1	6.2	18.1	6.6
10	28	18.3	16.7 - 17.2	17.5	6.4	18.5	6.8
11	29	19	17.3 - 17.9	17.9	6.6	18.9	7
12	30	19.7	18 - 18.5	18.3	6.8	19.3	7.2
12.5	31	20.3	18.6 - 19.1	18.7	6.9	19.7	7.3
13.5	32	21	19.2 - 19.8	19.1	7	20.1	7.4
1	33	21.7	19.9 - 20.5	19.5	7.1	20.5	7.5
2	34	22.3	20.6 - 21.1	19.9	7.3	20.9	7.7
2.5	35	23	21.2 - 21.8	20.3	7.4	21.3	7.8
3.5	36	23.7	21.9 - 22.5	20.7	7.6	21.7	8
4	37	24.3	22.6 - 23.1	21.1	7.7	22.1	8.1
5	38	25	23.2 - 23.8	21.5	7.9	22.5	8.3
6	39	25.7	23.9 - 24.5	21.9	8.1	22.9	8.5
6.5	40	26.3	24.6 - 25.1	22.3	8.3	23.3	8.7

Width 8		Anti Varus		DAFO		Shoe Size (UK)	Shoe Size (EU)
Width Circ. (cm)	Joint Circ. (cm)	Width Circ. (cm)	Joint Circ. (cm)	Width Circ. (cm)	Joint Circ. (cm)		
14.7	5.6	-	-	-	-	1	16
15.1	5.8	-	-	-	-	1.5	17
15.5	5.9	-	-	-	-	2.5	18
15.9	6	-	-	-	-	3	19
16.3	6.1	15.1	5.6	17.6	6.2	4	20
16.7	6.2	15.5	5.8	18	6.3	5	21
17.1	6.3	15.9	5.9	18.4	6.4	5.5	22
17.5	6.4	16.4	6.1	18.8	6.5	6	23
17.9	6.5	16.8	6.2	19.2	6.6	7	24
18.3	6.6	17.2	6.4	19.6	6.8	8	25
18.7	6.8	17.6	6.5	20	6.9	8.5	26
19.1	7	18	6.7	20.3	7	9.5	27
19.5	7.2	18.4	6.9	20.7	7.1	10	28
19.9	7.4	18.9	7	21.1	7.3	11	29
20.3	7.6	19.3	7.2	21.5	7.4	12	30
20.7	7.7	19.7	7.3	21.9	7.5	12.5	31
21.1	7.8	20.1	7.5	22.3	7.6	13.5	32
21.5	7.9	20.5	7.6	22.7	7.8	1	33
21.9	8.1	20.9	7.8	23.1	7.9	2	34
22.3	8.2	21.3	7.9	23.5	8	2.5	35
22.7	8.4	21.8	8.1	23.9	8.1	3.5	36
23.1	8.5	22.2	8.2	24.3	8.2	4	37
23.5	8.7	22.6	8.4	24.6	8.4	5	38
23.9	8.9	23	8.6	25	8.5	6	39
24.3	9.1	23.4	8.7	25.4	8.6	6.5	40