

## Aspen Vista® Universal Collar



### Features and Benefits

- Innovative height adjustment technology
- Adjustable front depth using a unique dial adjustment mechanism
- The Vista collar is six collars in one, allowing for improved fit
- Reducing storage and inventory costs
- 100% cotton - latex free

### Indications for Use

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- R/A & O/A cervical spine
- Herniated cervical disc
- Post-operative
- Motor Neurone Disease

Details	Size	Product Code
Aspen Vista Collar	Universal	984000
Aspen Vista Collar & Pad	Universal	984002

Details	Product Code
Replacement Pads (set)	984020
Back Panel Replacement Pads - Standard	984026
Back Panel Replacement Pads - Large	984028
Back Panel - Standard	984006
Back Panel - Large	984008

1

2

3

4

5

6

7

8

9