

## ToeOFF® 2.0 Splint



### Features and Benefits

- The ToeOFF 2.0 is the updated version of the ToeOFF Original
- Lightweight, carbon reinforced rehabilitation orthosis that absorbs force on heel strike
- Worn at the anterior side of the leg, bridging the ankle and preventing plantar flexion
- SoftKIT 2.0 comprises padding shaped to precisely fit the orthosis and comes supplied with two tibia pads
- ComfortKIT 2.0 comprises thicker padding, designed for the more sensitive skin, fits the orthosis precisely

### Indications for Use

- Mild to severe footdrop
- Mild to moderate ankle instability
- Impaired knee control
- Limb proprioception deficit
- Mild to moderate spasticity
- Partial foot amputation

**Note:** For ToeOFF Original please contact Customer Services.

Size	Footplate Dimensions			ToeOFF 2.0 Splint		Order ComfortKIT 2.0 and SoftKIT 2.0 Separately			
	Length (cm)	Height (cm)	UK Shoe Size	Product Code Left	Product Code Right	Size	Product Code SoftKIT 2.0	Product Code ComfortKIT 2.0	Replacement Straps
X-Small	21	36	1 - 3	289201010	289202010	X-Small	287500010	287510010	287600010
Small	23	38	3 - 4	289201011	289202011	Small	287500011	287510011	287600011
Medium	24.5	40.5	5 - 9	289201012	289202012	Medium	287500012	287510012	287600012
Large	27	43	10 - 11	289201013	289202013	Large	287500013	287510013	287600013
X-Large	28.5	45	12 - 13	289201014	289202014	X-Large	287500014	287510014	287600014

Note: The ToeOFF 2.0 when purchased individually does not come with padding. Please order this separately.

## BlueROCKER™ 2.0



### Features and Benefits

- The BlueROCKER 2.0 is the updated version of BlueROCKER Original
- BlueROCKER 2.0 is identical in shape to the ToeOFF 2.0 but offers more rigid orthotic control
- Developed primarily for bilateral patients and those with more involved pathologies
- The extra stability will improve both balance and posture and give the wearer greater security, especially individuals with weak quadricep muscles
- Usually the preferred orthosis to be used in conjunction with a socket and toe filler for management of partial foot amputations
- SoftKIT 2.0 comprises padding shaped to precisely fit the orthosis and comes supplied with two tibia pads
- ComfortKIT 2.0 comprises thicker padding, designed for the more sensitive skin, fits the orthosis precisely

### Indications for Use

- Footdrop
- Weakness or impairment in multiple leg muscle groups
- Impaired balance
- Limb proprioception
- Partial foot amputation

**Note:** For BlueROCKER Original please contact Customer Services.

Size	Footplate Dimensions		BlueROCKER 2.0		Order ComfortKIT 2.0 and SoftKIT 2.0 Separately			
	Length (cm)	Height (cm)	Product Code Left	Product Code Right	Size	Product Code SoftKIT 2.0	Product Code ComfortKIT 2.0	Replacement Straps
X-Small	21	36	289401010	289402010	X-Small	287500010	287510010	287600010
Small	23	38	289401011	289402011	Small	287500011	287510011	287600011
Medium	24.5	40.5	289401012	289402012	Medium	287500012	287510012	287600012
Large	27	43	289401013	289402013	Large	287500013	287510013	287600013
X-Large	28.5	45	289401014	289402014	X-Large	287500014	287510014	287600014

Note: The BlueROCKER 2.0 when purchased individually does not come with padding. Please order this separately.