

Calzamedi Footwear Style Options

Black Leather - Lace

Code: 4200-0101



Black Leather - Hook and loop fastening

Code: 4201-0101



Blue Nubuck

Hook and loop fastening 4203-0202

Lace 4202-0202



Pink Nubuck

Hook and loop fastening 4203-0203


Lace 4202-0203



Measuring for width:



Measuring for length:



Shoe Size		Inner Shoe Length	Foot Length	Width Circ.	4 Width	Width Circ.	6 Width	Width Circ.	8 Width	Width Circ.	Anti Varus Width	Width Circ.	DAFO Width
UK	EU												
1	16	10.3	9.2 - 9.7	12.7	4.8	13.7	5.2	14.7	5.6				
1.5	17	11	9.8 - 10.3	13.1	5	14.1	5.4	15.1	5.8				
2.5	18	11.7	10.4 - 10.9	13.5	5.1	14.5	5.5	15.5	5.9				
3	19	12.3	11 - 11.5	13.9	5.2	14.9	5.6	15.9	6				
4	20	13	11.6 - 12.1	14.3	5.3	15.3	5.7	16.3	6.1	15.1	5.6	17.6	6.2
5	21	13.7	12.2 - 12.8	14.7	5.4	15.7	5.8	16.7	6.2	15.5	5.8	18	6.3
5.5	22	14.3	12.9 - 13.3	15.1	5.5	16.1	5.9	17.1	6.3	15.9	5.9	18.4	6.4
6	23	15	13.4 - 14	15.5	5.6	16.5	6	17.5	6.4	16.4	6.1	18.8	6.5
7	24	15.7	14.1 - 14.7	15.9	5.7	16.9	6.1	17.9	6.5	16.8	6.2	19.2	6.6
8	25	16.3	14.8 - 15.2	16.3	5.8	17.3	6.2	18.3	6.6	17.2	6.4	19.6	6.8
8.5	26	17	15.3 - 15.9	16.7	6	17.7	6.4	18.7	6.8	17.6	6.5	20	6.9
9.5	27	17.7	16 - 16.6	17.1	6.2	18.1	6.6	19.1	7	18	6.7	20.3	7
10	28	18.3	16.7 - 17.2	17.5	6.4	18.5	6.8	19.5	7.2	18.4	6.9	20.7	7.1
11	29	19	17.2 - 17.9	17.9	6.6	18.9	7	19.9	7.4	18.9	7	21.1	7.3
12	30	19.7	18 - 18.5	18.3	6.8	19.3	7.2	20.3	7.6	19.3	7.2	21.5	7.4
12.5	31	20.3	18.6 - 19.1	18.7	6.9	19.7	7.3	20.7	7.7	19.7	7.3	21.9	7.5
13.5	32	21	19.2 - 19.8	19.1	7	20.1	7.4	21.1	7.8	20.1	7.5	22.3	7.6
1	33	21.7	19.9 - 20.5	19.5	7.1	20.5	7.5	21.5	7.9	20.5	7.6	22.7	7.8
2	34	22.3	20.6 - 21.1	19.9	7.3	20.9	7.7	21.9	8.1	20.9	7.8	23.1	7.9
2.5	35	23	21.2 - 21.8	20.3	7.4	21.3	7.8	22.3	8.2	21.3	7.9	23.5	8
3.5	36	23.7	21.9 - 22.5	20.7	7.6	21.7	8	22.7	8.4	21.8	8.1	23.9	8.1
4	37	24.3	22.6 - 23.1	21.1	7.7	22.1	8.1	23.1	8.5	22.2	8.2	24.3	8.2
5	38	25	23.2 - 23.8	21.5	7.9	22.5	8.3	23.5	8.7	22.6	8.4	24.6	8.4
6	39	25.7	23.9 - 24.5	21.9	8.1	22.9	8.5	23.9	8.9	23	8.6	25	8.5
6.5	40	26.3	24.6 - 25.1	22.3	8.3	23.3	8.7	24.3	9.1	23.4	8.7	25.4	8.6

Ordering Example (for Black Leather - Lace)

Select Size (EU) 32

Select Style (from top of page) 4200-0101

Select Width 4

= Ordering Code 32-4200-0101-4