

Comfort Night Splint



Features and Benefits

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and achilles tendon
- Soft lycra lined foam interior is soft, washable and breathable to eliminate heat build-up
- Lycra lined foam padded straps maintain comfort across the flexor tendons
- Lightweight, low profile shell is sturdy and breathable
- The orthoses helps to prevent pain and further injury by limiting contractions of the soft tissues at night
- Universal fit - left/right
- Adjustable flexion straps provide a gentle stretching of the plantar fascia and achilles tendon

Indications for Use

- Plantar fasciitis
- Achilles tendonitis
- Plantar flexion contractures

Size	UK Shoe Size	Product Code
Small	Up to 6	UFO/S/E
Medium	6 - 9	UFO/M/E
Large	9 +	UFO/L/E

Leeder Multi-Use Boot



Features and Benefits

- Utilises Flex-E-Core technology for adjustment without the need for tools or heat

Indications for Use

- Plantar fasciitis
- Achilles tendonitis
- Night splinting
- Foot drop

Size	Max. Foot Length (cm)	Brace Height (cm)	Product Code
Paediatric	15	25	LMUP
Adult	25	35	LMUA

Podus Boot



Features and Benefits

- A lightweight, quality plastic shell with poly/pile liner that helps treat and prevent lower extremity disorders associated with trauma or immobility
- An adjustable toe extension helps protect toes from bedding and the anti-rotation bar may be repositioned to assist with internal or external rotation
- Removable transfer pad on the dorsal side for limited patient ambulating
- Fits right or left foot
- Latex free

Indications for Use:

- Pressure neurosis
- Ankle contractures
- Post-op hip abduction/rotation

Size	Men's UK Shoe Size	Ladies' UK Shoe Size	Product Code
Regular	Up to 10	Up to 11	POD-R
X-Large	10.5 +	11.5 +	POD-XL

Replacement liners code: POD-Liners. Replacement liners available in pack of 3.