

Spinal Plus



Features and Benefits


- An ultralight and non-intrusive spinal brace for the treatment of osteoporosis
- Shoulder straps with pads that can be easily positioned by hook and loop fastening
 - Adjustable internal structure in aluminium alloy (radiotransparent)
 - Dorsal-abdominal stays with low-friction strap
 - Self-directing buckles

Indications for Use

- Lower back pain
- Osteoporosis
- Hyperkyphosis
- Juvenile osteochondrosis

Size	Distance T1 - L5 (cm)	Product Code
X-Small	35 - 41	MT1049XS
Small	41 - 47	MT1049S
Medium	47 - 57	MT1049M
Large	54 - 61	MT1049L

Thoraco Lumbar Support



Features and Benefits

- Highly ventilated elastic back panel with two removable shaped stays for increased support
- Stretch elastic cross over shoulder straps to maintain posture
- Touch and close fastening for ease of use

Indications for Use

- Osteoporosis - kyphosis reduction
- Thoraco lumbar muscular injuries

Size	Waist Circ. (cm)	Product Code
Small	61 - 70	TLS/S
Medium	71 - 81	TLS/M
Large	82 - 93	TLS/L
X-Large	94 - 104	TLS/XL
XX-Large	105 - 115	TLS/XXL
XXX-Large	115 - 127	TLS/XXXL

Universal Cross Pectoral Brace



Features and Benefits

A comfortable, cosmetic hyperextension brace which can easily be donned and doffed.

- Adjustable, flexible pectoral pads with silicone hinges fit comfortably below clavicles and does not impinge on the patient's sternum or choke when sitting, thus providing increased comfort and compliance
- ATL pectoral pads provide rotation and linear adjustment
- Hinged pubic pad assembly adjusts to the patient's contour and can be quickly manipulated for a precise fit
- Universal closure strap and pad
- One long strap which can be cut to the patient's circumference and fit with a Y-tab over the end
- Anterior closure allows for easy donning and doffing

Indications for Use

- Osteoporosis
- Kyphosis
- Thoracic stabilisation
- Thoracic extension
- Spinal alignment

Size	Description	Product Code
Universal	Universal Cross Pectoral Brace	CROSSP