

Diaped Duosoft Therapeutic



Features and Benefits

- Constructed with a multi-density material combination of plastazote and cellular urethane, these advanced insoles offer superb resilience, cushioning and shock absorption
- The insoles are designed to fit into both traditional footwear and 'additional depth' orthopaedic shoes
- The addition of an increased metatarsal support in the anatomical design is designed to facilitate off-loading of the metatarsal head area

Indications for Use

- For individuals with diabetes, arthritis, or other conditions that result in sensitive feet
- Metatarsalgia
- Pre-ulcerated sites
- Plantar fasciitis
- Foot deformities
- Morton's neuralgia
- Forefoot accommodation
- Skin breakdown

UK Shoe Size	EU Shoe Size	Product Code
2 - 4.5	34 - 37	SB2335
5 - 6.5	36 - 38	SB2336
7 - 8.5	39 - 42	SB2337
9 - 10.5	43 - 45	SB2338
11 - 13	46 - 48	SB2339

Slim Tech



Features and Benefits

- The lightweight and easy fitting Slim Tech with its 3° low arch allows for heel contact and a good comfortable fit
- The orthotic is for use in ladies' dress shoes and sports boots
- Slip is reduced by the design of the hook heel

Indications for Use

- Arch pain
- Knee pain
- Shin splints
- Tendon pain
- Back pain
- Postural related symptoms

UK Shoe Size	Product Code	UK Shoe Size	Product Code
3 - 4	SC3070	7 - 8	SC3072
5 - 6	SC3071	9 - 10	SC3073

Trio Slimfit



Features and Benefits

- The recessed medial heel post within the Trio Slimfit design offers control over pronation without risking irritation of the skin within the tighter fitting fashion shoe
- The forefoot contouring provides a stabilising effect upon the forefront joints and reduces cramping pain

Indications for Use

- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Size	UK Shoe Size	Product Code Medium Density
X-Small	3 - 4.5	TSF-XS
Small	5 - 6.5	TSF-S
Medium	7 - 8.5	TSF-M
Large	9 - 10.5	TSF-L