

## Equiflex



### Features and Benefits

- By elevating the heel the Equiflex will restore ankle joint extension reducing the demand for compensation
- Lateral forefoot is supported at the level of the metatarsal heads to facilitate this and accommodate any additional postings that may be required
- Supplied with 0°, 3° and 5° rear posts and 6mm elevators

### Indications for Use

- High arches
- Excessive pronation
- Forefoot pain
- Heel valgus

Size	UK Shoe Size	Product Code Medium Density	Product Code Firm Density
X-Small	3.5 - 4.5	EQM-XS	EQF-XS
Small	5 - 6	EQM-S	EQF-S
Medium	6.5 - 7.5	EQM-M	EQF-M
Large	8 - 9	EQM-L	EQF-L
X-Large	9.5 - 10.5	EQM-XL	EQF-XL
XX-Large	11 - 12	EQM-XXL	EQF-XXL

## Trio Easy Fit



### Features and Benefits

- The Easy Fit is intrinsically a 3/4 length medially posted orthotic with an additional posting option for triplanar control of pronation featuring a lateral cutaway design
- This model maintains the essential features for pronation control, including an optional biplanar rearfoot post
- A medial forefoot post can be used to replace the first metatarsal cushion if required
- Comfort is achieved without compromising natural movement and comfort

### Indications for Use

- Excessive pronation
- Forefoot pain
- Knee pain
- Heel valgus

Size	UK Shoe Size	Product Code Medium Density
Kids	1 - 2.5	TEF-K
X-Small	3 - 4.5	TEF-XS
Small	5 - 6.5	TEF-S
Medium	7 - 8.5	TEF-M
Large	9 - 10.5	TEF-L
X-Large	11 - 12.5	TEF-XL

## Trio Heel Pain



### Features and Benefits

- The Trio Heel Pain orthosis is contoured to support the navicular bone of the arch, directly reducing the cause of plantar fascia strain
- Higher loading forces under the heel can also lead to inflammation of the skin and supporting structures of the foot
- Features a large cushioning protection for the softer structure under the heel bones

### Indications for Use

- Plantar fasciitis
- Plantar heel pain
- Calf strain
- Anterior medial knee pain
- Hip instability

Size	UK Shoe Size	Product Code Medium Density
Small	5 - 6.5	THP-S
Medium	7 - 8.5	THP-M
Large	9 - 10.5	THP-L