

Escort Hip Brace



Features and Benefits

- A lightweight patient friendly brace using the latest in hip joint technology
- The joint system abducts the leg only during flexion allowing a normal gait
- The pelvic section is designed to fit just below the iliac crest, eliminating impingement on the rib cage whilst seated
- The thigh section is a two-piece medial/lateral design that allows for maximum adjustability and positive purchase on the leg

Indications for Use

- Functional and effective support in the treatment of various hip pathologies

Size	Waist Size (cm)	Thigh Circ. (cm)	Product Code Left	Product Code Right
Small	66 - 81	43 - 53	HIP-ES-SL	HIP-ES-SR
Medium	81 - 97	51 - 61	HIP-ES-ML	HIP-ES-MR
Large	97 - 112	58 - 69	HIP-ES-LL	HIP-ES-LR
X-Large	112 - 127	66 - 76	HIP-ES-XLL	HIP-ES-XLR

Arthritic Knee Sleeve



Features and Benefits

- Simple knee sleeve
- Pull on design, elasticated for easy use of donning
- Provides compression and warmth for painful joints
- Posterior shaping eradicates bunching in popliteal area

Indications for Use

- Arthritis

Size	Knee Circ. (cm)	Product Code
Small	23 - 26	803/AKS
Medium	26 - 29	805/AKS
Large	29 - 32	807/AKS
X-Large	32 - 35	808/AKS
XX-Large	35 - 39	809/AKS

Standard Knee Brace



Features and Benefits

- The front opening design allows easy application for the arthritic or elderly
- Metal side hinges provide medial/lateral stability while allowing flexion
- Opposing fastenings for upper and lower touch and close straps counter rotation and allow correct fitting
- Removable hinges
- Universal left/right application

Indications for Use

- Medio lateral stability
- Osteoarthritis
- Rheumatoid arthritis

Material: Airprene

Size	Knee Circ. (cm)	Product Code
X-Small	24 - 28	SKB/CX/BE/XS
Small	28 - 34	SKB/CX/BE/S
Medium	34 - 39	SKB/CX/BE/M
Large	39 - 44	SKB/CX/BE/L
X-Large	44 - 49	SKB/CX/BE/XL