

## Philadelphia Universal Cervical Collar w/Trachea



### Features and Benefits

- Featuring a two-piece design anatomically contoured to provide superior stability and comfort with large trachea opening for easy access
- Constructed from durable plastic, laminated to soft hypo-allergenic foam with contact closure

### Indications for Use

- Immobilisation for both extrication and rehabilitation
- Post-operative discectomy
- Spinal fusion

Height measurement taken from posterior portion of mandible to the clavicle directly below.

Size	Height (cm)	Neck Circ. at Adam's Apple (cm)	Product Code
Universal	5 - 7	20.3 - 28	19-4100

## Ambu Perfit Collar



### Features and Benefits

- Features one-piece flat design with flip out chin piece and trachea opening for easy access
- Constructed from plastic, laminated to foam with contact closure

### Indications for Use

- Immobilisation for both extrication and rehabilitation

Size	Colour	Product Code
Infant	Purple	19-9001
Pedi	Aqua	19-9002
Universal Mini Adjustable	No specific colour	19-9008

## Ultra Cervical Collar



### Features and Benefits

- Features anatomically contoured foam with soft cotton stockinette covering for maximum patient comfort
- Constructed from medium firm density foam with contact closure

### Indications for Use

- Cervical facet syndrome
- Discogenic pain
- Degenerated disc disease
- Muscle sprains or strains of cervical region

Size	Height (cm)	Neck Circ. at Adam's Apple (cm)	Product Code
XX-Small	5.1	25.4 - 30.5	15-5000
X-Small	5.1	35.6 - 40.6	15-5001

