

Basic Foam Collar



Features and Benefits

- Symmetrical design provides comfortable support, suitable for night use
- Touch and close fastening

Material: Medium density foam covered with stockinette.

Indications for Use

- Spondylitis
- Spondylosis
- Rheumatoid arthritis
- Osteoarthritis

Size	Neck Circ. (cm)	Front Depth (cm)	Product Code
X-Small	33 - 37	7	BFC/XS
Small	38 - 42	8	BFC/S
Medium	43 - 46	9	BFC/M
Large	47 - 51	10	BFC/L
X-Large	52 - 56	11	BFC/XL

*Made to measure service available, please provide neck circumference and front depth

Semi Rigid Collar



Features and Benefits

- Provides improved stability with shaped chin support, suitable for night use
- Touch and close fastening

Material: High density PVC foam laminate covered with stockinette.

Indications for Use

- Spondylitis
- Spondylosis
- Rheumatoid arthritis
- Osteoarthritis

Size	Neck Circ. (cm)	Front Depth (cm)	Product Code
Small	38 - 42	6.5	SRC/S
Medium	43 - 46	7	SRC/M
Large	47 - 51	8	SRC/L

*Made to measure service available, please provide neck circumference and front depth

Shaped Foam Collar



Features and Benefits

- Symmetrical design provides comfortable support, suitable for night use
- Increased shaping to reflect underlying anatomy and provide an improved fit
- Touch close fastening

Material: Medium density foam covered with stockinette.

Indications for Use

- Spondylitis
- Spondylosis
- Rheumatoid arthritis
- Osteoarthritis

Size	Neck Circ. (cm)	Front Depth (cm)	Product Code
Small	38 - 42	6.5	SFC/S
Medium	43 - 46	7	SFC/M
Large	47 - 51	8	SFC/L

*Made to measure service available, please provide neck circumference and front depth