

Foot-Up

1



Foot-Up shown with footwear

Features and Benefits

- Foot-Up is a lightweight ankle orthosis that offers dynamic support for drop foot or similar complaints
- The revolutionary orthosis provides visible improvement in gait by supporting the foot the moment it is raised
- Lightweight and difficult to see under clothes
- Constructed from a padded breathable material it is comfortable to wear with or without footwear (shoeless wrap required for use without footwear)

2

Indications for Use

- Peroneal neuropathies
- Drop foot
- Dysfunctional gait

3



Foot-Up shown with shoeless wrap

Foot-Up			
Size	Ankle Circ. (cm)	Product Code Black	Product Code Beige
1 - Medium	18 - 21	07810-1	07811-1
2 - Large	22 - 26	07810-2	07811-2
3 - X Large	27 - 33	07810-3	07811-3

Shoeless Wrap*			
Size	Foot Circ. (cm)	Product Code Black	Product Code Beige
1 - Small	< 20	27810-1	27811-1
2 - Medium	20 - 24	27810-2	27811-2
3 - Large	24 - 28	27810-3	27811-3
4 - X Large	> 28	27810-4	27811-4

*Shoeless Wrap comprises of foot wrap only and must be ordered with complete Foot-Up kit unless the patient already has one.

Foot-Up Plastic Insert		
Size	Product Code Black	Product Code Beige
All	17810	17811

Foot-Up and Shoeless Wrap Bundle*	
Product Code Black	Product Code Beige
07810-XY	07811-XY

*Where X = Size of Foot-Up & Y = Size of Shoeless Wrap, for example to order a Black size 1 (Medium) Foot-Up with a size 2 (Medium) Shoeless Wrap, the code would be 07810-12.

4

5

Universal PRAFO 650 Kodel

6



Features and Benefits

The PRAFO® (Pressure Relief Ankle Foot Orthosis) is a totally adjustable, custom-fitted AFO that can help manage many of the ankle/foot anomalies patients present. Its custom-contoured aluminium heel connector bar helps to control dorsi-plantar flexion by allowing infinite, measurable adjustments, whilst its ergonomically designed foot plate allows the liner to be changed easily without tools. Its integral, securely fastened vibram walking base provides a shoe-type surface for ambulatory patients. Design ensures that the patient is protected around the vulnerable areas of the foot and heel and will remain so during ambulation. Universal liners are also available to order - use code PRAFO-LINER.

The calf and foot segments can be custom-modified to accommodate patient physiology and the anti-rotation bar is made of malleable aluminium and adjusts to control unwanted extremity rotation.

- Designed for recumbent or ambulant protection of the heel and foot
- Zero pressure or shear force at the heel
- Universal sizing to reduce stock holding.
- Resists plantarflexion contractures developing
- Custom modifications are possible
- One version fits both left and right
- Calibrated adjustment of plantar/dorsiflexion
- Integral walking base - designed for ambulation
- Washable and replaceable kodel liner system
- Adjustable foot, toe and calf sections
- Aluminium heel connector for strength
- Rotation bar resists internal/external rotation when used in bed
- Suitable for day and night use

7

8



9

Size	Product Code
Universal	PRFO
Universal Liner	PRAFO-LINER

Indications for Use:

- Pressure ulcers
- Can be used to resist the development of contracture