

## Deviation Hand Orthosis



### Features and Benefits

- Terrycloth
- Headliner
- Comes with one liner
- Paediatric sizing
- Ambidextrous
- Progressively change wrist and finger positions and amount of wrist deviation by adjusting the internal frame and locking the frame's swivel with included hex wrench
- Correct or control wrist deviation motion even when changes in the position of the wrist and hand are made
- Unique below the wrist swivel ability that allows for lateral and medial movements
- Accommodates for ulnar and radial deviation

### Indications for Use

- Great for arthritic or neuromuscular deformities or for flaccid, weak or spastic hand conditions

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Paediatric Small	6.35 - 8.89	16.51	DH101-PS
Paediatric Medium	7.62 - 11.43	17.78	DH101-PM
Paediatric Large	11.43 - 15.24	22.86	DH101-PL

## Hand-Wrist Finger Orthosis



\*Shown in standard Terrycloth with finger separator attached



Available to order in Comfyprene purple colour, please contact Customer Services

### Features and Benefits

- Terrycloth
- Headliner
- Comes with one cover
- Ambidextrous
- Paediatric sizing
- Side wings prevent ulnar and radial deviation
- Arm cuffs easily adjust to width of arm
- Excellent resting splint to prevent trauma to joints and increases ROM
- Multi-layered for comfort
- Adjustable reinforced spine allows for easy adjustments without tools
- Easy to clean
- Keeps thumb in functional position without stretching thenar eminence
- Side wings prevent ulnar/radial deviation
- Washable and absorbent covers help prevent skin breakdown
- Soft straps accommodate fluctuations in edema
- Arm, wrist, fingers independently adjustable
- Resists strong flexion synergy

### Indications for Use

- This popular splint provides support and positioning for weak or deformed hands

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Paediatric Small	6.35 - 8.89	17.78	H101-PS
Paediatric Medium	7.62 - 11.43	20.32	H101-PM
Paediatric Large	11.43 - 15.24	24.13	H101-PL

\*Please note, the standard model supplied is Terrycloth.

**Note:** Supplied with finger separators