

Comfyrene Contour Hand Orthosis



Features and Benefits

- The Comfyrene® Contour Hand Orthosis has unique side flaps that hold the wrist in place allowing for full hand contouring
- C-bar and thumb support with included strap
- 4 straps come as standard
- Unique butterfly designed side straps hold wrist and forearm in place
- Contours the whole hand, wrist and forearm
- Instead of removable terrycloth covers, the Comfyrene splint inserts are covered with soft neoprene and headliner material, which creates a slim profile for everyday comfort
- Splints are adjustable and re-adjustable using the Comfy™ patented reinforced spine
- Straps may be positioned anywhere on the splint for optimal fit

Indications for Use

- Use after sports injuries
- Removal of casts
- Sports medicine
- Arthritic patient's support of the weakened extremity

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code Left	Product Code Right
Adult	17.8 +	27.9	HCO-101-CP-LB-LT	HCO-101-CP-LB-RT

Comfy Hand-Wrist Finger Orthosis



Features and Benefits

- Provides support and positioning for the weak or deformed wrist, hand or fingers
- Serves as an excellent resting splint to prevent trauma to joints and positions to increase ROM
- Side wings prevent ulnar and radial deviation
- Arm cuffs easily adjust to width of arm
- Excellent resting splint to prevent trauma to joints and increases ROM
- Terrycloth
- Headliner
- Ambidextrous
- Paediatric sizing

Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Carpal tunnel syndrome

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Small	15.2 - 17.8	27.9	H101-AS
Medium	17.8 - 20.3	31.8	H101-AM
Large	20.3 +	34.3	H101-AL

Add finger separators - FS1

Paediatric sizes also available - see page 86

Comfy Air Hand Orthosis



Features and Benefits

- The Comfy Air Hand Orthosis is an exceptionally useful splint for the severely contracted hand
- It has an adjustable air bladder that allows gradual extension of the MP, DIP and PIP joints, as well as an adjustable wrist for customised fit
- Ingeniously designed by a therapist, one simply inserts the un-fitted splint with the bladder under the fingers and then pumps air to inflate the bladder and increase finger extension. (Headliner material only)
- Ambidextrous

Indications for Use

- Ideal for the severely contracted hand
- Inflatable air bladder allows for gradual MP, DIP and PIP joint extension
- Adjustable wrist
- Provides low pressure stretch therapy
- Hand air bladder can be used separately without full orthosis as needed
- Can fit into even the most tightly clenched hand

Size	Wrist to Tip of Middle Finger (cm)	Total Length (cm)	Product Code
Adult	17.8 - 20.3	31.8	HA101-A

Short Wrist Brace



Features and Benefits

Standard single compression one-piece wrist brace provides good conformity around the wrist and hand.

Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Carpal tunnel syndrome

Material: Neoprene

		Product Code - Neoprene Beige	
Size	Wrist Circ. (cm)	Left	Right
X-Small	10 - 12.5	SWB/E/XS/L	SWB/E/XS/R
Small	12.5 - 15	SWB/E/S/L	SWB/E/S/R
Medium	15 - 17.5	SWB/E/M/L	SWB/E/M/R
Large	17.5 - 20	SWB/E/L/L	SWB/E/L/R
X-Large	20 - 22.5	SWB/E/XL/L	SWB/E/XL/R

Short Wrist Brace - Elastic



Features and Benefits

Standard single compression one-piece wrist brace provides good conformity around wrist and hand.

Indications for Use

- Flaccid drop wrist
- Repetitive strain injury
- Rheumatoid arthritis
- Carpal tunnel syndrome

Material: Elastic

		Product Code - Elastic Beige	
Size	Wrist Circ. (cm)	Left	Right
X-Small	10 - 12.5	SWB/XSL	SWB/XSR
Small	12.5 - 15	SWB/SL	SWB/SR
Medium	15 - 17.5	SWB/ML	SWB/MR
Large	17.5 - 20	SWB/LL	SWB/LR
X-Large	20 - 22.5	SWB/XLL	SWB/XLR