



## Mini TLC™ Hip Abduction Brace



### Features and Benefits

- Prefabricated
- Material: MPE
- Pelvic band adjustable for hip development and pelvic circumference with quick release buckle closure
- Thigh cuff bivalve design adjustable to thigh proximal and distal measurements
- Breathable, washable liner absorbs perspiration and reduces heat build up
- Universal hip joint, TLC hip joint with infinite abduction/adduction & flexion/extension settings
- Available in 6 sizes

### Indications for Use

- Arthroplasty
- Hip dislocation or potential high risk dislocation
- Hip dislocation - anterior or posterior
- Hip management/immobilisation
- Hip revision
- Hip surgery
- Post-operative
- Pre-operative
- Stabilise, align & reinforce hip

\*Available to order

Size	Waist Circ. (cm)	Proximal Thigh Circ. (cm)	Distal Thigh Circ. (cm)	Product Code
XX-Small	40.6 - 45.7	29.2 - 34.3	20.3 - 25.4	*TLCMINI2S
X-Small	45.7 - 50.8	31.8 - 36.8	22.9 - 27.9	*TLCMINIXS
Small	50.8 - 55.9	34.3 - 39.4	25.4 - 30.5	*TLCMINISM
Medium	55.9 - 61	36.8 - 41.9	27.9 - 33	*TLCMINIMD
Large	61 - 66	39.4 - 44.4	30.5 - 35.6	*TLCMINILG
X-Large	66 - 71.1	41.9 - 47	33 - 38.1	*TLCMINIXL

### Adjustable Abduction Bar

Size	Length of Bar (cm)
Juvenile	34.3
Paediatric	30.5