

Paediatric Coloured Neoprene Wrist/Thumb Brace



Features and Benefits

- Manufactured in a range of coloured neoprene
- Contrasting colourful bar casing and hook and loop fastening
- 15cm length
- Sold in single units
- Colours available: black, royal blue, navy, beige, red, day-glo orange

Indications for Use

- Healing fractures
- Healing ligament injuries

Size	Wrist Circ. (cm)	Product Code Left	Product Code Right
Small	Up to 11.5	HL-PNT-1	HR-PNT-1
Medium	11.5 - 14.5	HL-PNT-2	HR-PNT-2
Large	14.5 - 16	HL-PNT-3	HR-PNT-3

Pavlik Harness



Features and Benefits

- Strong cotton webbing with comfort padding
- Hook and loop fastening with colour coded tabs
- Soft cotton bootees
- Effective in the treatment of congenital dislocation of the hips

Indications for Use

- Indicated for babies with suspected congenital dislocation of the hip (CDH)

Size	Chest Circ. (cm)	Product Code
X-Small	38 - 44.5	PHXS
Small	44.5 - 53	PHS
Medium	53 - 60	PHM
Large	61 - 67	PHL
X-Large	68 - 77	PHXL

Paediatric Hip Abduction Splint - Standard & Adjustable



Features and Benefits - Standard Model

- Maintains concentric reduction of hips in patients with hip dysplasia
- Replaces spica casts and harnesses for certain developmental dysplasia patients
- One-piece, non-adjustable polypropylene orthosis holds infant in 90° hip flexion and 60° abduction
- Open design facilitates nappy changes
- Closed-cell foam lining does not retain moisture

Indications for Use - Standard Model

- Indicated for patients with hip dysplasia

Features and Benefits - Adjustable Model

- Lightweight, paediatric orthosis offers advantage of abduction adjustment over the standard paediatric hip splint model described above
- Sturdy polypropylene frame set in 90° hip flexion - abduction easily adjustable 40° through 180° by aluminium slotted screw assembly
- Fits infants from 3 months (6-12 lbs) to 18 months (over 30 lbs).
- More total contact than Pavlik harness
- Easy application, removal facilitates nappy changes
- Closed-cell foam does not retain moisture

Indications for Use - Adjustable Model

- Indicated for patients with hip dysplasia
- X-Large size effective for certain cerebral palsy patients

Size	Average Age	Weight Guide in (lbs)	Product Code Standard Model	Product Code Adjustable Model
X-Small	Up to 3 mths	6 - 12	B/3899	B/3906
Small	3 - 6 mths	13 - 16	B/3900	B/3907
Medium	6 - 12 mths	17 - 21	B/3901	B/3908
Large	1 year	22 - 30	B/3902	B/3909
X-Large	18 mths	30+	B/3903	B/3910