


Hyperextension Brace



Features and Benefits

- Lightweight aluminium construction
- Padded pelvic band - can be locked or free to adjust as patient changes from standing to sitting, avoiding brace shift
- Simple height and width adjustment of sternal component
- Three step auto-lock system allows for easier access
- Lightweight and hard wearing polymer hinge with no metal springs guarantees the right amount of movement between the frames and pelvic band


Indications for Use

- Immobilisation following inter-vertebral disc surgery
- Stable vertebral compression fractures of the lower thoracic and lumbar spine

Size	Pelvis Circ. (cm)	Height of Brace (cm)*	Product Code
X-Small	50 - 60	30 - 37	MP100XS
Small	60 - 75	34 - 42	MP100S
Medium Short	75 - 90	34 - 41	MP100MS
Medium	75 - 90	36 - 46	MP100M
Medium Long	75 - 90	40 - 52	MP100ML
Large Short	90 - 105	36 - 46	MP100LS
Large	90 - 105	40 - 52	MP100L
X-Large Short	105 - 115	39 - 49	MP100XLS
X-Large	105 - 115	42 - 54	MP100XL

*Measured from the proximal side of the body of the sternum to the upper edge of the pubis (pubic tubercle)

SOMI Brace



Features and Benefits

Fully adjustable cervical system providing maximum support whilst immobilising the head, includes headband for rigidity when chin plate is removed for washing, eating or dressing.

Indications for Use

- Post trauma cervical spine instability
- Rheumatoid arthritis control of extension/flexion - side flexion and rotation

Size	Product Height (cm)	Product Code
Small	Up to 144	SOMI/SML
Regular	Over 145	SOMI/REG