


## Escort Hip Brace



### Features and Benefits

- A lightweight patient friendly brace using the latest in hip joint technology
- The joint system abducts the leg only during flexion allowing a normal gait
- The pelvic section is designed to fit just below the iliac crest, eliminating impingement on the rib cage whilst seated
- The thigh section is a two-piece medial/lateral design that allows for maximum adjustability and positive purchase on the leg

### Indications for Use

- Functional and effective support in the treatment of various hip pathologies

Size	Waist Size (cm)	Thigh Circ. (cm)	Product Code Left	Product Code Right
Small	66 - 81	43 - 53	HIP-ES-SL	HIP-ES-SR
Medium	81 - 97	51 - 61	HIP-ES-ML	HIP-ES-MR
Large	97 - 112	58 - 69	HIP-ES-LL	HIP-ES-LR
X-Large	112 - 127	66 - 76	HIP-ES-XLL	HIP-ES-XLR

## Arthritic Knee Sleeve



### Features and Benefits

- Simple knee sleeve
- Pull on design, elasticated for easy use of donning
- Provides compression and warmth for painful joints
- Posterior shaping eradicates bunching in popliteal area

### Indications for Use

- Arthritis

Size	Knee Circ. (cm)	Product Code
Small	23 - 26	803/AKS
Medium	26 - 29	805/AKS
Large	29 - 32	807/AKS
X-Large	32 - 35	808/AKS
XX-Large	35 - 39	809/AKS

## Standard Knee Brace



### Features and Benefits

- The front opening design allows easy application for the arthritic or elderly
- Metal side hinges provide medial/lateral stability while allowing flexion
- Opposing fastenings for upper and lower touch and close straps counter rotation and allow correct fitting
- Removable hinges
- Universal left/right application

### Indications for Use

- Medio lateral stability
- Osteoarthritis
- Rheumatoid arthritis

Material: Airprene

Size	Knee Circ. (cm)	Product Code
X-Small	24 - 28	SKB/CX/BE/XS
Small	28 - 34	SKB/CX/BE/S
Medium	34 - 39	SKB/CX/BE/M
Large	39 - 44	SKB/CX/BE/L
X-Large	44 - 49	SKB/CX/BE/XL

## Long Knee Brace



### Features and Benefits

- Longer version of the standard knee brace (40cm) provides greater medial/lateral support and support for the varus/valgus knee
- Universal left/right application

### Indications for Use

- Medio lateral stability and patella control osteoarthritis
- Osteoarthritis
- Rheumatoid arthritis

Material: Airprene

Size	Knee Circ. (cm)	Product Code
Small	27 - 35.5	LKBB/CX/BE/S
Medium	33 - 40.5	LKBB/CX/BE/M
Large	38 - 45.5	LKBB/CX/BE/L
X-Large	43 - 55.5	LKBB/CX/BE/XL

## Stabilised Knee Sleeve



### Features and Benefits

- Incorporates medial and lateral flexible spiral stays for enhanced medial/lateral stability
- Provides compression and support
- Touch and close straps for improved fit
- Universal left/right application

### Indications for Use

- Medio lateral stability and patella control knee instability
- Post-op rehabilitation
- Strains/sprains
- Bursitis
- Tendonitis osteoarthritis
- Rheumatoid arthritis
- Osteoarthritis

Material: Airprene

Size	Knee Circ. (cm)	Product Code
X-Small	22 - 28	SKS/CX/BE/XS
Small	28 - 34	SKS/CX/BE/S
Medium	34 - 39	SKS/CX/BE/M
Large	39 - 44	SKS/CX/BE/L
X-Large	44 - 49	SKS/CX/BE/XL
XX-Large	49 - 56	SKS/CX/BE/XXL

## Deluxe Knee Immobiliser



### Features and Benefits

- Lightweight three panel design with easy application due to removable side panels
- Malleable steels for contouring to the patient
- Universal left/right application

### Indications for Use

- Post trauma immobilisation
- Post-op immobilisation

Length (cm)	Product Code
35	DKI/350
40	DKI/400
45	DKI/450
50	DKI/500
55	DKI/550
60	DKI/600